RR0210=00000_K, RR0310=00000_K

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1.1 Danger during Preparation for Use

⚠ CAUTION

Risk of injury due to incorrect environmental conditions for the electric strap rewinder.

The electric strap rewinder (option) as part of the standing frame is a MEDICAL ELECTRIC DEVICE according to EN 60601-1-2. To exclude the risk of injuries when using the electric strap rewinder, the following must be observed:

- The product may only be installed by qualified personnel.
- During preparation for use, special safety precautions apply in regards to electromagnetic compatibility (EMC).
- Note that the device may only be installed and put into operation subject to the EMC information contained in this accompanying document.

2 Tables

2.1 Electromagnetic emission

Guidance and manufacturer's declaration – electromagnetic emission		
The STRAP REWINDER FOR STANDING FRAMES is intended for use in the electromagnetic environment specified below. The customer or the user of the STRAP REWINDER FOR STANDING FRAMES should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The STRAP REWINDER FOR STANDING FRAMES uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The STRAP REWINDER FOR STANDING FRAMES is suitable for use in all establishments, including domestic establishments and those directly connected to the public lowvoltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

2.2 Electromagnetic immunity

Guidance and manufacturer's declaration – electromagnetic immunity			
The STRAP REWINDER FOR STANDING FRAMES is intended for use in the electromagnetic environment specified below. The customer or the user of the STRAP REWINDER FOR STANDING FRAMES should assure that it is used in such an environment			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (>95 % dip in U_T) for 5 sec	< 5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the STRAP REWINDER FOR STANDING FRAMES requires continued operation during power mains interruptions, it is recommended that the STRAP REWINDER FOR STANDING FRAMES be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a. c. mains voltage prior to application of the test level.			

2.3 Electromagnetic immunity

Guidance and manufacturer's declaration – electromagnetic immunity			
The STRAP REWINDER FOR STANDING FRAMES is intended for use in the electromagnetic environment specified below. The customer or the user of the STRAP REWINDER FOR STANDING FRAMES should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the STRAP REWINDER FOR STANDING FRAMES, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1,17\sqrt{P}$ $d=1,17\sqrt{P}$ 80 MHz to 800 MHz $d=2,34\sqrt{P}$ 800 MHz to 2,5 GHz where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). ^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the STRAP REWINDER FOR STANDING FRAMES is used exceeds the applicable RF compliance level above, the STRAP REWINDER FOR STANDING FRAMES should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the STRAP REWINDER FOR STANDING FRAMES.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

2.4 Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT

Recommended separation distances between portable and mobile RF communications equipment and the STRAP REWINDER FOR STANDING FRAME.

The STRAP REWINDER FOR STANDING FRAMES is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the STRAP REWINDER FOR STANDING FRAMES can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the STRAP REWINDER FOR STANDING FRAMES as recommended below, according to the maximum output power of the communications equipment

Rated maximum output of transmitter W	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d=1,17\sqrt{P}$	80 MHz to 800 MHz $d=1,17\sqrt{P}$	800 MHz to 2,5 GHz $d=2,34\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,69	3,69	7,38
100	11,67	11,67	23,33

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Otto Bock South Africa (Pty) Ltd
Building 3 Thornhill Office Park · 94 Bekker Road
Midrand · Johannesburg · South Africa
T +27 (0) 11 312 1255
info-southafrica@ottobock.co.za
www.ottobock.co.za

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Ottobock has a certified Quality Management System in accordance with ISO 13485.

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Otto Bock Korea HealthCare Inc.
Otto Bock South East Asia Co., Ltd.
Otto Bock HealthCare GmbH
Otto Bock Manufacturing Königsee GmbH
Otto Bock Mobility Solutions GmbH