


ottobock.

446B1=1.5-7; 446B1=2-7; 481C08=SK079; 481
S47=SK001 / 481S47=SK002 / 481C17=SK020

DE	Gebrauchsanweisung	4
EN	Instructions for Use	7
IT	Istruzioni per l'uso	10
ES	Instrucciones de uso	13
SV	Bruksanvisning	16

INFORMATION

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
INFORMATION

Date of the last update: 2013-09-03

- Please read this document carefully.
- Follow the safety instructions.

Explanation of Symbols

 Warnings regarding possible risks of accident or injury.

 Information regarding operation. Information for service personnel.

1 Description

1.1 Intended Use

The lap belt with snap buckle may be used **exclusively** to support a correct seating and pelvis position of the user in the wheelchair or seating shell.

1.2 Safety Instructions

INFORMATION

This lap belt does not replace the safety belt when transporting patients in a vehicle.

Risk of burns in the proximity of fire. Falling embers can make supports, rests, plates, shells, pads, covers or cushions begin to smoulder and thus cause fire.

Keep away from all sources of ignition, especially lit cigarettes.

1.3 Function

The purpose of the lap belt is to provide posture support by positioning the pelvis of the patient in the seating system. Lap belts should be used

with all seating systems. They are essential when torso or head positioning systems are employed.

The belts are available in two widths - 35 mm and 50 mm. They are length-adjustable in order to assure correct and secure positioning of the patient.

The 481S47=SK001 / 481S47=SK002 / 481C17=SK020 and 446B1=1.5-7 / 446B1=2-7 lap belts (Figure 1 or similar) can be mounted to the wheelchair frame with the help of D-rings or loops.

The 481C08=SK079 lap belt (Figure 2) can only be mounted to the wheelchair frame with loops.

Lap belt components (Figure 3):

1. Belt with snap buckle

2. Glide buckle

only 481S47=SK001 / 481S47=SK002 / 481C17=SK020; 446B1=1.5-7 / 446B1=2-7:

3. D-ring

4. Cap screw (length varies depending on the wheelchair type)

5. M6 washer

6. M6 nut

7. Round-head self-tapping screw

2 Handling

2.1 Mounting on the Wheelchair System

Attaching the D-ring to the lap belt: Thread the belt (Figure 4, item 1) through the glide buckle (Figure 4, item 2), through the D-ring (Figure 4, item 3) and back through the glide buckle (Figure 4).

Ensure that the loose end of the belt faces out (Figure 4, item 4) so that subsequent adjustments can be made easily.

Mounting the D-ring on the wheelchair frame: See Figure 3, item 4 - 6 for the correct positioning of the screw, nut and washer. Use existing holes or screw connections for mounting where possible.

Round-head self-tapping screws are included with the belts for use on seating shells etc. (Figure 3, item 7).

Mounting the lap belt using loops: The lap belt can be mounted by passing belt loops around the frame components and attaching them (Figure 4, item 5).

2.2 Adjustment

Open the lap belt by pressing on the top and bottom halves of the snap buckle simultaneously and position the patient in the wheelchair.

The lap belt is adjusted via the glide buckles. Both halves of the belt should be adjusted to the same length so that the closure is centred.

The lap belt should be tight enough to provide sufficient support for the pelvis, but not so tight as to be uncomfortable. The adjustment angle between the lap belt and seat frame should be approximately 45° - 60°.

3 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used properly and maintained according to the instructions for use.

The manufacturer is not responsible for damages caused by components and spare parts not approved by the manufacturer. Repairs must be carried out exclusively by authorised dealers or by the manufacturer.

4 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock Mobility Solutions GmbH