Cloud 476C00=FFC*

DE Gebrauchsanweisung	3
EN Instructions for Use	8
FR Instructions d'utilisation	13
IT Istruzioni per l'uso	18
ES Instrucciones de uso	23
PT Manual de utilização	28
NL Gebruiksaanwijzing	33
SV Bruksanvisning	38
DA Brugsanvisning	43
FI Käyttöohje	48
PL Instrukcja użytkowania	53
CS Návod k použití	58
TR Kullanma talimatı	63

INFORMATION

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INFORMATION

Date of the last update: 2014-02-03

- Please read this document carefully.
- Follow the safety instructions.

Explanation of Symbols

CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Information regarding operation. Information for service personnel.

1 Description**1.1 Intended Use**

The **Cloud** wheelchair seat cushion is to be used **exclusively** and in dependence of the degree of severity for decubitus prophylaxis and decubitus therapy for patients with paraplegia, hemiplegia, neurological deficits, rheumatic nervous and muscle disorders, as well as for geriatric patients, and is to be used exclusively in a wheelchair.

For hygienic reasons, re-use of the seat cushion or use by another patient is not recommended by the manufacturer.

1.2 Field of Application

The seat cushion can be used by patients with a body weight of up to 125 kg / 275 lbs without restrictions.

CAUTION

Risk of redness of the skin/decubitus when exceeding the maximum permissible patient weight. When exceeding the maximum permissible patient weight, we highly recommend that the user test the cushion for a few hours under medical supervision, during which time it should be determined that no reddening of the skin occurs, before using the cushion on a permanent basis.

1.3 Contraindication

Relative contraindications: Should redness of the skin occur, immediately stop using the seat cushion and consult the treating physician or therapist.

1.4 Safety Instructions

The **Cloud** wheelchair seat cushion has been manufactured in accordance with state-of-the-art technology and is safe in operation. The safety of the seat cushion has been confirmed by CE certification and the declaration of conformity. The development, construction and production of the product completely meet the requirements of Medical Device Directive 93/42/EEC.

1.4.1 General Safety Instructions

INFORMATION

Seat cushions may only be selected and adapted by appropriately trained and qualified personnel.

⚠ CAUTION

Risk of tipping. Inserting a cushion into a wheelchair may change the wheelchair's centre of gravity, which can cause the wheelchair in some cases to tip over backwards and may result in injuries. It is therefore important to check whether anti-tippers or adjustable axle supports, such as those used for lower limb amputees, should be mounted to the wheelchair to improve stability.

⚠ CAUTION

Risk of pressure sores. All Ottobock cushions have been designed to reduce pressure sores. A cushion alone, however, cannot completely eliminate sitting pressure or possible pressure sores. Ottobock cushions are not intended to replace thorough skin care, good nutrition, hygiene and regular pressure relief.

⚠ CAUTION

Risk of redness of the skin. Before using a cushion on a permanent basis, we highly recommend that the user test it for a few hours under medical supervision, during which time it should be determined that no reddening of the skin occurs. Examine the skin for redness on a regular basis. Redness of the skin is a clinical indicator for tissue damage. Should redness of the skin occur, immediately stop using the cushion and consult a physician or therapist.

NOTICE

Risk of damaging the surface of the cushion. Do not allow sharp objects to come into contact with the cushion and do not expose the cushion to open flames.

NOTICE

Risk of damage due to impermissible temperature. Do not expose the cushion to temperatures of more than 55 °C / 131 °F for a longer period of time.

If the cushion was exposed to temperatures of less than 0 °C / 32 °F for a longer period of time, touch the cells to test how they feel. If they feel hard, knead the cells with your hands until they are fluffy again. Place the cushion on the wheelchair.

1.5 Function

The **Cloud** cushion is tailored to the special needs of active users and offers a maximum of protection. The use of Floam polymer gel represents a new way of decubitus prophylaxis. The Floam filling in the cushion provides for effective pressure distribution to lower-risk areas. The side cells of Floam provide additional support.

Cloud consists of:

- Cover (1)
- Foam side cells (2)
- Individual Foam cells (3)
- Foam base (4)

Sizes

The seat area is divided into 16 cells which are anatomically arranged on a foam base:

- Blue Foam cell: small
- Black Foam cell: medium
- Red Foam cell: large

Cel-Pak

The inner cover (Cel-Pak) serves to hold the side cells of Floam and the individual Floam cells in position and tightly encloses the foam base.

After removal of the foam base from the bottom side of the cushion, qualified personnel can exchange or replace any individual Floam cell, if necessary.

1.6 Accessories

The seat cushion can be combined with an incontinence cover (article number 477C00=ST410/..411/..412/..413).

2 Handling

2.1 Positioning

INFORMATION

Pressure on the cushions can only be optimally distributed, if the footrest of the wheelchair has been adjusted in such a way that the thigh of the user is in full contact with the cushion surface.

Always make sure that the cushion is correctly positioned in the wheelchair and that the cover is correctly re-mounted after it has been removed.

Place the cushion on the wheelchair. The **“back”** label attached to the bottom side of the cushion should point towards the backrest. The square hook and loop patch sewn to the bottom side of the cushion makes the cushion easy to fix in place.

2.2 Function Test

Daily:

Loosen up the **Cloud** cushion by pushing the polymer gel of the Floam cells from the sides towards the centre of the cushion. Do not take the cushion out of the cover.

Since the Floam's own weight does not it to flow, this measure makes sure that a maximum of Floam is underneath the user.

Monthly:

For inspection of the Floam cells you can remove the cushion cover.

Remove the foam basis on the bottom side of the cushion and then remove the cells through the slots (5).

Inspect the individual Floam cells. If necessary, they can be replaced by qualified personnel.

2.3 Assembly

For cleaning it may be necessary to remove the cushion parts from the cover. To do so, open the zipper of the cover and remove the cushion parts.

For reassembly, please proceed in reverse order paying attention to the following:

- Insert the foam parts/Floam pads through the slots of the inner cover (Cel-Pak) (2/3).

INFORMATION

When inserting the Floam cells into the Cel-Pak (inner cover) make sure that they match the location according to their colour.

- The foam base covers the inner cover (Cel-Pak) (4).
- Pull the outer over the inner cover (Cel-Pak) so that the foam base faces down, the Floam pad faces up and the “back” labels show in the same direction.

2.4 Cleaning and Care

NOTICE

Risk of damage due to improper cleaning. Do not use any bleaching agents. Do not tumble-dry the foam parts.

INFORMATION

Should it be necessary to disinfect the cushion (cover/foam parts), add a standard hygiene detergent during washing (e.g. Sagrotan laundry hygiene detergent).

In case of heavy incontinence, an additional incontinence cover should be used. Wash the cover separately at 40 °C / 104 °F using a mild, environment-friendly detergent.

Cover

For separate washing, the zipper of the cover can be opened and the cover removed.

1. Turn the cover inside out for cleaning.
 2. The cover can be washed at 65 °C / 149 °F using a mild, environment-friendly detergent. To prevent excessive wear we recommend, however, to wash the cover at 40 °C / 104 °F (gentle wash).
 3. Allow to air dry. When using an automatic clothes drier, select the gentle cycle at low temperature.
- Replace the cover in case of advanced wear and tear.

Foam parts

Hand-wash the foam parts at 40 °C / 104 °F using a mild, environment-friendly detergent. Allow to air dry.

Floam cells

Should it be necessary in rare cases to clean the Floam cells, then please observe the following: The sheath of the Floam cells is made of double-layer polyurethane foil. For cleaning, remove it from the Cel-Pak cover. Rinse it with lukewarm water, or use a disinfectant or detergent that are not based on alcohol.

Dry the cells with a towel or let them air dry at room temperature.

3 Technical Data

3.1 Available Sizes

The seat cushion will be adapted to the correct size in accordance with the order:

Cloud seat cushion	Dimensions seat width x seat depth
[] 476C00=FFC14x16	35x40
[] 476C00=FFC15x15	38x38
[] 476C00=FFC15x16	38x40
[] 476C00=FFC16x16	40x40
[] 476C00=FFC16x16,5	40x42
[] 476C00=FFC16x18	40x46
[] 476C00=FFC16,5x16	42x40

Cloud seat cushion	Dimensions seat width x seat depth
[] 476C00=FFC16,5x16,5	42x42
[] 476C00=FFC17x18	43x46
[] 476C00=FFC17,5x16,5	44x42
[] 476C00=FFC18x16	46x40
[] 476C00=FFC18x16,5	46x42
[] 476C00=FFC18x18	46x46
[] 476C00=FFC19x18	48x46

3.2 Basic Data

	Seat cushion thickness*	Weight*	Permissible patient weight	Reference size*
Cloud	11 cm	3.2 kg (with 42 x 40 cm/ 16.5 x 16.5 inches)	125 kg / 275 lbs	42 x 40 cm (16.5 x 16 inches)
	Structure of the seat cushions			
	Foam: Lycra, Denier, Pennyla, Velcro, urethane Foam cells: Floam (mineral oil, Kraton, acrylic) Cover: decor fabric, black			

* ready-to-use state

4 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

The manufacturer is not responsible for damages caused by components and spare parts not approved by the manufacturer. Repairs must be carried out exclusively by authorised dealers or by the manufacturer.

5 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock Korea HealthCare Inc.
South East Asia Co., Ltd.
Otto Bock HealthCare GmbH