

A200

Power wheelchair

Max. load capacity: 100 kg/220 lbs

Article number: 490E46=00000_C



Standard Equipment

- VR2 control unit
- Control panel holder, standard
- 6 km/h
- 28 Ah (C5) AGM batteries
- Battery charger
- 8" PU caster wheels
- Standard seat
- Seat cushion
- Armrests
- Footrests
- Silver metallic frame colour


Order no.:	

<input type="checkbox"/> Provide estimate only <input type="checkbox"/> Order	
<input type="checkbox"/> Attachment for custom fabrication enclosed	
Billing address/customer number	
Delivery address	
Sales person	Name Date Signature

Basic data (voluntary information)

Gender _____ Indication _____

Age _____ Paying party _____



Mandatory information*

Body height in cm _____

Body weight in kg _____

1 Seat depth _____

2 Trunk depth _____

3 Armrest height _____

4 Armpit height _____

5 Shoulder height _____

6 Head height _____

7 Forearm length to the middle of the palm _____

8 Trunk width _____

9 Lower leg length _____

10 Seat width _____

* This information will help to support user-oriented settings and product development. Ottobock guarantees the confidential and purely internal use of these data.

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please use the spare parts catalogue for spare parts orders.

Information




EC 01 Standard seat



EC 03 Contour seat



ED 20 Mechanical back angle adjustment 30°



EP 65 Lap belt with buckle



EB 02 Individual footrests with aluminium footplates



EB 03 Single-panel footrest



EB 20 Mechanically elevating footrest



EB 80 Lateral leg pad



EF 07 Curb climbing assist, to overcome curbs with a maximum height of 10 cm

Standard seat																			
EC 18	<input type="checkbox"/>	Standard seat Junior (SD 34 ... 40 cm, SW 34 ... 40 cm)	Add. charge																
EC 01	<input type="checkbox"/>	Standard seat, small (SD 38 ... 46 cm, SW 38 ... 42 cm)	-																
EC 02	<input type="checkbox"/>	Standard seat, large (SD 42 ... 50 cm, SW 43 ... 48 cm)	-																
Settings for the Standard seat																			
Back angle	<input type="checkbox"/>	-9/1/11/21°	<input type="checkbox"/>	0/10/20/30°															
Not available with mechanical back angle adjustment!																			
Seat depth SD in cm	<input type="checkbox"/>	34	<input type="checkbox"/>	36	<input type="checkbox"/>	38	<input type="checkbox"/>	40	<input type="checkbox"/>	42	<input type="checkbox"/>	44	<input type="checkbox"/>	46	<input type="checkbox"/>	48	<input type="checkbox"/>	50	
Seat width SW cm																		
Lower leg length LL cm (25 ... 47 cm) without seat cushion																		
Seat height SH in cm	<input type="checkbox"/>	43	<input type="checkbox"/>	48	(Note: without seat cushion, dependent on lower leg length, seat height less lower leg length equals ground clearance, a minimum of 3 cm is recommended. Seat height 43 cm not available with standard seat, large)														
Contour seat																			
EC 03	<input type="checkbox"/>	Contour seat, small (SD 36 ... 44 cm, SW 38 ... 42 cm)	Add. charge																
EC 04	<input type="checkbox"/>	Contour seat, large (SD 40 ... 48 cm, SW 43 ... 48 cm)	Add. charge																
Settings for the contour seat																			
Seat depth SD in cm	<input type="checkbox"/>	36	<input type="checkbox"/>	38	<input type="checkbox"/>	40	<input type="checkbox"/>	42	<input type="checkbox"/>	44	<input type="checkbox"/>	46	<input type="checkbox"/>	48					
Seat width SW cm																		
Lower leg length LL cm (28 ... 54 cm)																		
Seat height SH in cm	<input type="checkbox"/>	approx. 51	<input type="checkbox"/>	approx. 56	(Note: (with seat padding, dependent on LL) SH less LL equals ground clearance, a minimum of 3 cm is recommended) SH of 51 cm not available with contour seat, large!														
Seat cushion (only for Standard seat)																			
EC 06	<input type="checkbox"/>	Seat cushion, black cover (standard equipment)	-																
EC 10	<input type="checkbox"/>	Seat cushion with incontinence cover	Add. charge																
EC 07	<input type="checkbox"/>	Without seat cushion	Price reduction																
EC 11	<input type="checkbox"/>	Terra seat cushion (up to 140 kg/309 lbs)	Add. charge																
EC 12	<input type="checkbox"/>	Terra Aquos seat cushion (up to 140 kg/309 lbs)	Add. charge																
EC 13	<input type="checkbox"/>	Aerius seat cushion (up to 140 kg/309 lbs)	Add. charge																
EC 14	<input type="checkbox"/>	Aerius Active seat cushion (up to 150 kg/330 lbs)	Add. charge																
EC 17	<input type="checkbox"/>	Terra Flair seat cushion (up to 150 kg/330 lbs)	Add. charge																
Back adjustment																			
ED 20	<input type="checkbox"/>	Mechanical back angle adjustment	Add. charge																
Footrests																			
EB 01	<input type="checkbox"/>	Individual footrests, with plastic footplates (standard equipment) Not available with Standard seat Junior.	-																
EB 02	<input type="checkbox"/>	Individual footrests, with aluminium footplates (continuously angle-adjustable)	Add. charge																
EB 03	<input type="checkbox"/>	Footrest with continuous aluminium footplate (continuously angle-adjustable) Not available with Standard seat Junior.	Add. charge																
EB 07	<input type="checkbox"/>	Without footrest	Price reduction																
EB 20	<input type="checkbox"/>	Mechanically elevating footrests with length adjustment and aluminium footplates, angle adjustable (LL with standard seat from 27 cm, with contoured seat from 32 cm)	Add. charge																
Footrest accessories																			
EB 80	<input type="checkbox"/>	Lateral leg pad for mechanically elevating footrest (pair)	Add. charge																
EB 70	<input type="checkbox"/>	Fixation straps for aluminium footplates (included with mounting accessories)	Add. charge																

Armrest accessories			
EE 34	<input type="checkbox"/>	Channel forearm pad adapter, fixed, mounting side left	Add. charge
EE 24	<input type="checkbox"/>	Channel forearm pad adapter, fixed, mounting side right	Add. charge
EE 35	<input type="checkbox"/>	Channel forearm pad adapter, swivelling, mounting side left	Add. charge
EE 25	<input type="checkbox"/>	Channel forearm pad adapter, swivelling, mounting side right	Add. charge
EE 36	<input type="checkbox"/>	Arm strap, mounting side left	Add. charge
EE 26	<input type="checkbox"/>	Arm strap, mounting side right	Add. charge

For information on channel forearm pads and hand pads, please refer to our accessories catalogue.

Frame colour			
016	<input type="checkbox"/>	Silver metallic (standard equipment)	-
038	<input type="checkbox"/>	Sparkle red	Add. Charge
004	<input type="checkbox"/>	Ice blue	Add. Charge
001	<input type="checkbox"/>	Brilliant silver	Add. Charge
002	<input type="checkbox"/>	Anthracite metallic	Add. Charge
025	<input type="checkbox"/>	Stainless steel silver	Add. Charge
018	<input type="checkbox"/>	Signal red RAL 3001	Add. Charge
003	<input type="checkbox"/>	Jet black RAL 9005	Add. Charge
009	<input type="checkbox"/>	Vinho sparkle	Add. Charge
010	<input type="checkbox"/>	Marine blue	Add. Charge
021	<input type="checkbox"/>	Sparkle granny-smith	Add. Charge
005	<input type="checkbox"/>	Sparkle light orange	Add. Charge
020	<input type="checkbox"/>	Saturn blue metallic	Add. Charge
027	<input type="checkbox"/>	Shimano matte	Add. Charge
017	<input type="checkbox"/>	Traffic yellow RAL 1023	Add. Charge
012	<input type="checkbox"/>	Sky blue	Add. Charge
014	<input type="checkbox"/>	Candy red	Add. Charge
015	<input type="checkbox"/>	Clear reptile glossy	Add. Charge
013	<input type="checkbox"/>	Jet black RAL 9005 matte	Add. Charge
019	<input type="checkbox"/>	White pearl	Add. Charge
026	<input type="checkbox"/>	Sparkling iron effect light	Add. Charge
028	<input type="checkbox"/>	Rainbow	Add. Charge

Tyres			
EF 10	<input type="checkbox"/>	Puncture-proof tyres for caster wheels (standard equipment)	-
EG 06	<input type="checkbox"/>	Puncture-proof tyres for all wheels	Add. charge

Control panel holder			
EU 19	<input type="checkbox"/>	Standard control panel holder (standard equipment)	-
EU 18	<input type="checkbox"/>	Control panel holder, swing-away, height-adjustable	Add. charge
EU 20	<input type="checkbox"/>	Control panel holder, swing-away in parallel	Add. charge

Control panel holder mounting side			
	<input type="checkbox"/>	Right	
	<input type="checkbox"/>	Left	

Control unit accessories			
EU 15	<input type="checkbox"/>	Attendant control	Add. charge
EU 50	<input type="checkbox"/>	Handheld programming device (dealer version, not subject to discount)	Add. charge
EU 55	<input type="checkbox"/>	Programming software, includes cable (dealer version, not subject to discount, certificate required)	Add. charge

EU 23	<input type="checkbox"/>	Fork for tetraplegics, horizontal (adjustable for hand width 70 ... 90 mm)	Add. charge
EU 24	<input type="checkbox"/>	Stick S80 (special joystick handle 80 mm)	Add. charge
EU 25	<input type="checkbox"/>	Soft ball (black cover)	Add. charge
EU 26	<input type="checkbox"/>	Ball top, small (diameter 43 mm)	Add. charge
EU 27	<input type="checkbox"/>	Ball top, large (diameter 48 mm)	Add. charge
EU 28	<input type="checkbox"/>	Joystick top, flexible, including ball top, large	Add. charge
EU 35	<input type="checkbox"/>	Soft ball, coloured (diameter 30 mm, 4 pieces)	Add. charge

(Extension by 60 mm; control panel holder, downwards adjustable or swing-away, is recommended)

General accessories

EE 30	<input type="checkbox"/>	Adapter for headrest installation kit Please refer to our accessories catalogue for information about headrests.	Add. charge
EP 65	<input type="checkbox"/>	Lap belt with buckle	Add. charge
ES 25	<input type="checkbox"/>	Additional lights	Add. charge
EF 07	<input type="checkbox"/>	Curb climbing assist (not available in combination with single-panel footrest)	Add. charge
EZ 52	<input type="checkbox"/>	Tool kit with tyre pump	Add. charge
ES 20	<input type="checkbox"/>	Pannier bag, detachable	Add. charge
ES 22	<input type="checkbox"/>	Accessories bag (e.g. for tool kit)	Add. charge
ES 24	<input type="checkbox"/>	Pocket for mobile phone	Add. charge
ES 35	<input type="checkbox"/>	Tray, swing-away to the side (incl. mounting accessories; cut-out for control panel must be adapted)	Add. charge
EA 60	<input type="checkbox"/>	Crutch holder	Add. charge
ES 60	<input type="checkbox"/>	Vehicle transportation kit (required for AMF anchor point system), (only in combination with max 100 kg/220 lbs load capacity), (not available in combination with mechanical back angle adjustment or mechanically elevating footrests)	Add. charge

Special options

24817.001	<input type="checkbox"/>	Special colour two-layer coating	Add. charge
24817.000	<input type="checkbox"/>	Special colour single-layer coating	Add. charge
21187.002	<input type="checkbox"/>	Special dimensions for platform for seating shell	Add. charge
20465.016	<input type="checkbox"/>	Standard back tube extension/shortening	Add. charge
40405.002	<input type="checkbox"/>	Back upholstery for power wheelchairs with special length	Add. charge
20465.011	<input type="checkbox"/>	Contour back with standard seat bottom or contour seat bottom with standard back	Add. charge
36362.001	<input type="checkbox"/>	Avantgarde2 footplates, locked	Add. charge
25142.001	<input type="checkbox"/>	Amputation legrest for power wheelchairs	Add. charge
38068.000	<input type="checkbox"/>	A200 with enAble-40 control unit	Add. charge
25856.002	<input type="checkbox"/>	A200 base frame straight	Add. charge