

**Base price: EUR x,xxx**

All prices in EUR,  
valid from 1/1/2014 to 31/12/2014

Order no.:

\_\_\_\_\_

Provide estimate only     Order

Attachment for custom fabrication enclosed

Billing address/customer number

Delivery address

Sales person

Name

.....  
Date                      Signature

## Skippi plus

### Power wheelchair for children

Max. load capacity: 50 kg/110 lbs

**Article number: 490E55=00000\_K**




#### Standard Equipment

- enAble50 with the ability to control four electric options
- 6 km/h
- 28 Ah (C5) AGM batteries
- Battery charger
- 8" PU caster wheels
- Standard seat with back bag
- Seat cushion
- Telescoping side panels
- Individual telescoping footrests
- Wheel cover for drive wheels
- Drive wheels with aluminium rim, black
- Lap belt with buckle

#### Basic data (voluntary information)

Gender \_\_\_\_\_ Indication \_\_\_\_\_

Age \_\_\_\_\_ Paying party \_\_\_\_\_



#### Mandatory information\*

Body height in cm \_\_\_\_\_

Body weight in kg \_\_\_\_\_

**1** Seat depth \_\_\_\_\_

**2** Trunk depth \_\_\_\_\_

**3** Armrest height \_\_\_\_\_

**4** Armpit height \_\_\_\_\_

**5** Shoulder height \_\_\_\_\_

**6** Head height \_\_\_\_\_

**7** Forearm length to the middle of the palm \_\_\_\_\_

**8** Trunk width \_\_\_\_\_

**9** Lower leg length \_\_\_\_\_

**10** Seat width \_\_\_\_\_

\* This information will help to support user-oriented settings and product development. Ottobock guarantees the confidential and purely internal use of these data.

**The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please use the spare parts catalogue for spare parts orders.**

## Information



EE 37 Lateral pads



EC 50 Electric seat tilt



EB 02 Individual footrests




EU 15 Attendant control



EB 05 Single-panel footrest



EB 20 Mechanically elevating footrests



ES 63 Bumper bar

<b>Settings for standard seat</b>		(SD 28 ... 38 cm, SW 30 ... 38 cm)	
<b>Back angle</b>	<input type="checkbox"/> -9/1/11/21°	<input type="checkbox"/> 0/10/20/30°	
Not available with electric back angle adjustment!			
<b>Seat depth SD in cm</b>	<input type="checkbox"/> 28	<input type="checkbox"/> 30	<input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38
<b>Seat width SW</b>	..... cm (30 ... 38 cm, SW 24 ... 29 cm possible with removable lateral pads)		
<b>Lower leg length LL</b>	..... cm (15 ... 38 cm) <b>without seat cushion</b>		
<b>Seat options</b>			
<b>EC 51</b>	<input type="checkbox"/>	<b>Mechanical seat tilt</b> (height-adjustable, removable push handles required)	<b>Add. charge</b>
<b>EC 50</b>	<input type="checkbox"/>	<b>Electric seat tilt</b>	<b>Add. charge</b>
<b>Back adjustment</b>			
<b>ED 50</b>	<input type="checkbox"/>	<b>Electric back angle adjustment 25°</b> (not available in combination with seating shell adapter)	<b>Add. charge</b>
<b>Seating system accessories</b>			
<b>EC 57</b>	<input type="checkbox"/>	<b>Electric seat height adjustment 20 cm</b>	<b>Add. charge</b>
<b>EE 30</b>	<input type="checkbox"/>	<b>Mounting accessories for headrest attachment</b> Please refer to our accessories catalogue for information about headrests.	<b>Add. charge</b>
<b>EP 65</b>	<input type="checkbox"/>	<b>Lap belt with buckle</b> (standard equipment)	-
<b>EE 37</b>	<input type="checkbox"/>	<b>Removable lateral pad for side panels</b> (set, 3 cm thick) (not available in combination with seating shell adapter for Mygo seating unit)	<b>Add. charge</b>
<b>ED 10</b>	<input type="checkbox"/>	<b>Push handles, height-adjustable, removable</b>	<b>Add. charge</b>
<b>ED 11</b>	<input type="checkbox"/>	<b>Back extension</b> (35 cm – 45 cm)	<b>Add. charge</b>
<b>Footrests</b>			
<b>EB 02</b>	<input type="checkbox"/>	<b>Individual footrests, with aluminium footplates</b> (standard equipment, LL 23 ... 38 cm)	-
<b>EB 05</b>	<input type="checkbox"/>	<b>Single-panel footrest with aluminium footplate</b> (LL 23 ... 38 cm)	<b>Add. charge</b>
<b>EB 06</b>	<input type="checkbox"/>	<b>Footrest for short LL</b> (LL 15 ... 22 cm)	<b>Add. charge</b>
<b>EB 07</b>	<input type="checkbox"/>	<b>Without footrest</b>	<b>Price reduction</b>
<b>EB 16</b>	<input type="checkbox"/>	<b>Footrest for short lower leg length with single-panel aluminium footplate</b> (LL 15 ... 22 cm)	<b>Add. charge</b>
<b>EB 20</b>	<input type="checkbox"/>	<b>Mechanically elevating footrest</b> (LL 15 ... 33 cm)	<b>Add. charge</b>
<b>Frame colour</b>			
<b>017</b>	<input type="checkbox"/>	<b>Traffic yellow RAL 1023</b> (back upholstery blue, wheel covers red)	-
<b>018</b>	<input type="checkbox"/>	<b>Signal red RAL 3001</b> (back upholstery blue, wheel covers yellow)	-
<b>010</b>	<input type="checkbox"/>	<b>Marine blue</b> (back upholstery red, wheel covers yellow)	-
<b>016</b>	<input type="checkbox"/>	<b>Silver metallic</b> (back upholstery black, wheel covers black)	-
<b>Tyres</b>			
<b>EF 10</b>	<input type="checkbox"/>	<b>Puncture-proof tyres on the caster wheels</b> (standard equipment)	-
<b>EG 06</b>	<input type="checkbox"/>	<b>Puncture-proof tyres for all wheels</b>	<b>Add. charge</b>

Control unit/control panel		
EU 85	<input type="checkbox"/>	<b>Standard control panel</b> (standard equipment) -
EU 86	<input type="checkbox"/>	<b>Control panel with two external buttons</b> <b>Add. charge</b> Please indicate button type! <input type="checkbox"/> <b>Piko button</b> (diameter 50 mm, red/black) <input type="checkbox"/> <b>Piko button</b> (diameter 30 mm, red/black)

Control panel holder		
EU 19	<input type="checkbox"/>	<b>Standard control panel holder</b> (standard equipment) -
EU 20	<input type="checkbox"/>	<b>Control panel holder, swing-away in parallel</b> <b>Add. charge</b>

Control panel holder mounting side (not for mid-tray control)		
	<input type="checkbox"/>	<b>Right</b>
	<input type="checkbox"/>	<b>Left</b>

Mid-tray control		
EU 87	<input type="checkbox"/>	<b>Mid-tray control (MTC) electrically swivelling</b> (charging/programming box required) <b>Add. charge</b> Please indicate button type! <input type="checkbox"/> <b>Piko button</b> (diameter 50 mm, red/black) <input type="checkbox"/> <b>Piko button</b> (diameter 30 mm, red/black)

Tray-side mounting materials for mid-tray control (MTC)		
	<input type="checkbox"/>	<b>Without</b>
EU 871	<input type="checkbox"/>	<b>Tray with holder</b> <b>Add. charge</b>
EU 873	<input type="checkbox"/>	<b>Holder only</b> <b>Add. charge</b>

External buttons for mid-tray control (MTC)		
	<input type="checkbox"/>	<b>Without</b>
EU 874	<input type="checkbox"/>	<b>Piko button, diameter 30 mm, red</b> <b>Add. charge</b>
EU 875	<input type="checkbox"/>	<b>Piko button, diameter 30 mm, yellow</b> <b>Add. charge</b>
EU 876	<input type="checkbox"/>	<b>Piko button, diameter 30 mm, black</b> <b>Add. charge</b>
EU 877	<input type="checkbox"/>	<b>Piko button, diameter 50 mm, red</b> <b>Add. charge</b>
EU 878	<input type="checkbox"/>	<b>Piko button, diameter 50 mm, yellow</b> <b>Add. charge</b>
EU 879	<input type="checkbox"/>	<b>Piko button, diameter 50 mm, black</b> <b>Add. charge</b>

One button, "swivel", required; connection of max. 3 buttons, "swivel", "ON/OFF" and "mode" possible.

Charging/programming box for mid-tray control (MTC) (includes 1.5 m BUS cable)		
EU 56	<input type="checkbox"/>	<b>Yes</b> <b>Add. charge</b>

Special controls		
EU 07	<input type="checkbox"/>	<b>Push-button module</b> (control of max. 10 electric options without LCD) <b>Add. charge</b> (not available in combination with chin control)
EU 97	<input type="checkbox"/>	<b>Programmable LCD</b> , IR environment control, Bluetooth and LED control included <b>Add. charge</b>
EU 98	<input type="checkbox"/>	<b>Programmable LCD</b> , IR environment control and LED control included <b>Add. charge</b>
EY 10	<input type="checkbox"/>	<b>ASL interface</b> (incl. standard holder, also select code EU 75) (for USA only) <b>Add. charge</b>

The following special controls require an LCD:		
<b>EU 09</b>	<input type="checkbox"/>	<b>Chin control with mini joystick, electrically swing-away</b> <b>Add. charge</b> (standard control panel included; recommended for control with lips or tongue, the mini joystick is very sensitive and not suitable for strong control movements of the chin!) Please indicate: <input type="checkbox"/> <b>Soft ball with elastomer</b> <input type="checkbox"/> <b>Soft ball with spring, hard</b> <input type="checkbox"/> <b>Soft ball with spring, soft</b> Please indicate installation version: <input type="checkbox"/> <b>Chin swivel arm, left</b> <input type="checkbox"/> <b>Chin swivel arm, right</b>
	<input type="checkbox"/>	<b>Chin control with midi joystick, electrically swing-away</b> <b>Add. charge</b> (standard control panel included) Please indicate installation version: <input type="checkbox"/> <b>Chin swivel arm, left</b> <input type="checkbox"/> <b>Chin swivel arm, right</b>
<b>EU 11</b>	<input type="checkbox"/>	<b>Sip and puff control</b> (standard control panel included) <b>Custom fabrication</b>
<b>EU 40</b>	<input type="checkbox"/>	<b>Calibratable joystick, standard, includes multifunction button</b> <b>Add. charge</b> (standard control panel included)
<b>EU 45</b>	<input type="checkbox"/>	<b>Calibratable joystick incl. multifunction button for assembly in tray</b> <b>Add. charge</b> (standard control panel included)
Joystick top version: Please indicate.		
	<input type="checkbox"/>	<b>Control knob with elastomer</b>
	<input type="checkbox"/>	<b>Soft ball with elastomer</b>
	<input type="checkbox"/>	<b>Soft ball with spring, hard</b>
	<input type="checkbox"/>	<b>Soft ball with spring, soft</b>
<b>EU 59</b>	<input type="checkbox"/>	<b>Adapter cable for left/right mouse click</b> <b>Add. charge</b> please select: <input type="checkbox"/> <b>Piko button</b> (diameter 50 mm, red/black) <input type="checkbox"/> <b>Piko button</b> (diameter 30 mm, red/black)
The following options require the IR LCD display code EU 97 or EU 98:		
<b>EU 61</b>	<input type="checkbox"/>	<b>IR remote control hands-free telephone including mains adapter</b> <b>Add. charge</b>
<b>EU 63</b>	<input type="checkbox"/>	<b>"Dragger 32" mouse click input assistance</b> <b>Add. charge</b> (included for LCD with Bluetooth in combination with special controls)
<b>Special control accessories</b>		
<b>EU 68</b>	<input type="checkbox"/>	<b>Adapter for 1 – 6 buttons</b> (required for 1 – 6 button control) <b>Add. charge</b>
<b>EU 75</b>	<input type="checkbox"/>	<b>ON/OFF button for LCD module, includes adapter cable</b> <b>Add. charge</b> Not available in combination with calibratable joystick or chin control. Required for ASL interface. ON/OFF button required for LED control.
Piko button (diameter 50 mm)		
<b>EU 691</b>	<input type="checkbox"/>	<b>Red</b> <b>Add. charge</b>
<b>EU 693</b>	<input type="checkbox"/>	<b>Black</b> <b>Add. charge</b>
<b>EU 692</b>	<input type="checkbox"/>	<b>Yellow</b> <b>Add. charge</b>
<b>EU 697</b>	<input type="checkbox"/>	<b>Silver</b> <b>Add. charge</b>
Piko button (diameter 30 mm)		
<b>EU 701</b>	<input type="checkbox"/>	<b>Red</b> <b>Add. charge</b>
<b>EU 702</b>	<input type="checkbox"/>	<b>Yellow</b> <b>Add. charge</b>
<b>EU 703</b>	<input type="checkbox"/>	<b>Black</b> <b>Add. charge</b>
<b>EU 704</b>	<input type="checkbox"/>	<b>Silver</b> <b>Add. charge</b>

Activation of the drive-away lock		
	<input type="checkbox"/>	Yes
	<input type="checkbox"/>	No

Activation of the electronic steering lock		
	<input type="checkbox"/>	Yes
	<input type="checkbox"/>	No

LCD holder		
EU 31	<input type="checkbox"/>	Standard LCD holder (standard equipment for LCD) -
EU 30	<input type="checkbox"/>	LCD holder, swing-away, height- and angle-adjustable Add. charge
EU 32	<input type="checkbox"/>	LCD holder with gooseneck Add. charge
EY 34	<input type="checkbox"/>	LCD holder for tray (included) Add. charge

Control unit accessories		
EU 15	<input type="checkbox"/>	Attendant control Add. charge
EU 22	<input type="checkbox"/>	Protective bow for standard control panel Add. charge
EU 23	<input type="checkbox"/>	Fork for tetraplegics, horizontal (adjustable for hand width 70 ... 90 mm) Add. charge
EU 24	<input type="checkbox"/>	Stick S80 (special joystick handle 80 mm) Add. charge
EU 25	<input type="checkbox"/>	Soft ball (black cover) Add. charge
EU 26	<input type="checkbox"/>	Ball top, small (diameter 43 mm) Add. charge
EU 35	<input type="checkbox"/>	Soft ball, coloured (diameter 30 mm, 4 pieces) Add. charge
EU 50	<input type="checkbox"/>	Handheld programming device (dealer version, not subject to discount) Add. charge
EU 55	<input type="checkbox"/>	enAble programming software with USB adapter (not subject to discount, certificate required!) Add. charge
EU 681	<input type="checkbox"/>	2-fold adapter cable with Piko button (diameter 50 mm, black) Add. charge (free assignment of 2 electrical functions) Please indicate function: <input type="checkbox"/> Chin swivel arm <input type="checkbox"/> Seat tilt <input type="checkbox"/> Back angle adjustment <input type="checkbox"/> Seat height adjustment <input type="checkbox"/> Light <input type="checkbox"/> Warning flasher <input type="checkbox"/> Direction indicator, left <input type="checkbox"/> Direction indicator, right <input type="checkbox"/> Horn
EU 682	<input type="checkbox"/>	2-fold adapter cable with Piko button (diameter 30 mm, black) Add. charge (free assignment of 2 electrical functions) Please indicate function: <input type="checkbox"/> Chin swivel arm <input type="checkbox"/> Seat tilt <input type="checkbox"/> Back angle adjustment <input type="checkbox"/> Seat height adjustment <input type="checkbox"/> Light <input type="checkbox"/> Warning flasher <input type="checkbox"/> Direction indicator, left <input type="checkbox"/> Direction indicator, right <input type="checkbox"/> Horn

<b>General accessories</b>			
<b>ES 25</b>	<input type="checkbox"/>	<b>Additional lights</b>	<b>Add. charge</b>
<b>ES 05</b>	<input type="checkbox"/>	<b>Lighting according to German Motor Vehicle Safety Standards (StVZO)</b>	<b>Add. charge</b>
<b>ES 39</b>	<input type="checkbox"/>	<b>Manual horn</b>	<b>Add. charge</b>
<b>ES 35</b>	<input type="checkbox"/>	<b>Tray, swing-away to the side</b> (includes mounting accessories; cut-out for control panel must be adapted)	<b>Add. charge</b>
<b>ES 40</b>	<input type="checkbox"/>	<b>Pennant</b>	<b>Add. charge</b>
<b>ES 63</b>	<input type="checkbox"/>	<b>Bumper bar, removable</b> (only up to LL 32 cm)	<b>Add. charge</b>
<b>EZ 52</b>	<input type="checkbox"/>	<b>Tool kit with tyre pump</b>	<b>Add. charge</b>
<b>ES 60</b>	<input type="checkbox"/>	<b>Vehicle transportation kit</b> (required for AMF anchor point system) Note exclusion list, contact Ottobock if necessary	<b>Add. charge</b>