

## B500 advanced

### Power wheelchair E 504

Maximum load capacity: 140 kg/309 lbs  
**Article number: 490E32=00000\_C**



Figure shows add. charge options

Order no.: _____	
---------------------	--


<input type="checkbox"/> <b>Provide estimate only</b>		<input type="checkbox"/> <b>Order</b>	
<input type="checkbox"/> <b>Attachment for custom fabrication enclosed</b>			
Billing address/customer number			
Delivery address			
Sales person		Name	
		.....	
		Date	Signature

### Standard Equipment

- 6 km/h
- enAble50 impulse control with up to 4 electric adjustments
- 63 Ah (C5) gel batteries
- Battery charger
- Standard seat
- Seat cushion
- Armrests
- Footrests
- Control panel with LCD and buttons for lights
- Lap belt with buckle

### Basic data (voluntary information)\*

Gender \_\_\_\_\_ Indication \_\_\_\_\_  
 Age \_\_\_\_\_ Paying party \_\_\_\_\_



\* This information will help to support user-oriented settings and product development. Ottobock guarantees the confidential and purely internal use of these data.

### Mandatory information\*

Body height in cm \_\_\_\_\_  
 Body weight in kg \_\_\_\_\_  
 1 Seat depth \_\_\_\_\_  
 2 Trunk depth \_\_\_\_\_  
 3 Armrest height \_\_\_\_\_  
 4 Armpit height \_\_\_\_\_  
 5 Shoulder height \_\_\_\_\_  
 6 Head height \_\_\_\_\_  
 7 Forearm length to the middle of the palm \_\_\_\_\_  
 8 Trunk width \_\_\_\_\_  
 9 Lower leg length \_\_\_\_\_  
 10 Seat width \_\_\_\_\_

**The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please use the spare parts catalogue for spare parts orders.**

**Information**



EC 01/EC 02 Standard seat



EC 03/EC 04 Contour seat



EC 40 - EC 46 Recaro seating system



ED 50 Electric back angle adjustment 30°



EC 50 Electric seat tilt 20°



EC 54 Electric seat tilt 45°



EB 01 Individual footrests



EB 03 Single-panel footrest



EB 20 Mechanically elevating footrest



ES 50 Electric LED lighting




ET 04 Curb climbing assist



EG 05 Drive wheel suspension



EU 23 - EU 28 Joystick options



EA 62 Bumper



EB 80 Lateral leg pad



EP 51 Four-way chest strap



EU 07 Push-button module

Standard seat			
<b>EC 18</b>	<input type="checkbox"/>	<b>Standard seat, Junior</b> (SD 34 ... 40 cm, SW 34 ... 40 cm)	<b>Add. charge</b>
<b>EC 01</b>	<input type="checkbox"/>	<b>Standard seat, small</b> (SD 38 ... 46 cm, SW 38 ... 42 cm)	-
<b>EC 02</b>	<input type="checkbox"/>	<b>Standard seat, large</b> (SD 42 ... 50 cm, SW 43 ... 48 cm)	-
<b>EC 19</b>	<input type="checkbox"/>	<b>Standard seat, XL</b> (SD 42 ... 50 cm, SW 50 ... 56 cm)	<b>Add. charge</b>
Note: Standard seat, XL is not in combination with SH 45 cm possible. Standard seat, XL is not in combination with footrests possible: Optioncode: EB 01, EB 02, EB 20, EB 51.			
Settings for the standard seat			
<b>Back angle</b>		<input type="checkbox"/> -9/1/11/21° <input type="checkbox"/> 0/10/20/30°	
Note: Not available with mechanical and electric back adjustment. Back angle -9/1/11/21° is not possible in combination with SD 34 and 36 cm.			
<b>Seat inclination SI</b>		<input type="checkbox"/> 0° <input type="checkbox"/> 3° <input type="checkbox"/> 6°	
<b>EC 65</b>	<input type="checkbox"/>	<b>Seat inclination 20°</b> (only for standard seat, SH 45 and 50 cm)	<b>Add. charge</b>
Seat inclination of 6° and 20° not available with mechanical or electric seat tilt.			
<b>Seat depth SD in cm</b>		<input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50	
<b>Seat width SW</b>		..... cm	
<b>Lower leg length LL</b>		..... cm (25 ... 54 cm) <b>without seat cushion</b>	
<b>Seat height SH in cm</b>		<input type="checkbox"/> 45 <input type="checkbox"/> 50 <input type="checkbox"/> 55	
Note: In case of 3° seat inclination, without seat cushion, dependent on lower leg length; seat height less lower leg length equals ground clearance, a minimum of 6 cm is recommended.			
Contour seat			
<b>EC 03</b>	<input type="checkbox"/>	<b>Contour seat, small</b> (SD 36 ... 44 cm, SW 38 ... 42 cm)	<b>Add. charge</b>
<b>EC 04</b>	<input type="checkbox"/>	<b>Contour seat, large</b> (SD 40 ... 48 cm, SW 43 ... 48 cm)	<b>Add. charge</b>
<b>EC 21</b>	<input type="checkbox"/>	<b>Contour seat, XL</b> (SD 40 ... 48 cm, SW 50 ... 56 cm) (not in combination with seat height 52 cm possible)	<b>Add. charge</b>
Settings for the contour seat			
<b>Seat depth SD in cm</b>		<input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48	
<b>Seat width SW</b>		..... cm	
<b>Lower leg length LL</b>		..... cm (28 ... 54 cm)	
<b>Seat height SH in cm</b>		<input type="checkbox"/> approx. 52 <input type="checkbox"/> approx. 57	
(Notice: In case of 0° seat inclination, with seat padding, dependent on lower leg length; seat height less lower leg length equals ground clearance, a minimum of 6 cm is recommended.)			
Recaro seat			
<b>EC 40</b>	<input type="checkbox"/>	<b>Recaro n-Joy</b> (manually adjustable backrest)	<b>Add. charge</b>
<b>EC 41</b>	<input type="checkbox"/>	<b>Recaro n-Joy-Plus</b> (adaptable on site, manually adjustable backrest) <b>Seat depth SD 38 ... 56 cm</b>	<b>Add. charge</b>
Modular Recaro seating system			
<b>EC 42</b>	<input type="checkbox"/>	<b>Backrest LX</b> (manually adjustable)	<b>Add. charge</b>
<b>EC 43</b>	<input type="checkbox"/>	<b>Backrest LT</b> (lateral supports with different setting possibilities, manually adjustable)	<b>Add. charge</b>
<b>EC 44</b>	<input type="checkbox"/>	<b>Seat bottom F</b> (flat padded, seat bottom sewn to lateral guides)	<b>Add. charge</b>
<b>EC 45</b>	<input type="checkbox"/>	<b>Seat bottom X</b> (contoured lateral guides with seat extension)	<b>Add. charge</b>
<b>EC 46</b>	<input type="checkbox"/>	<b>Seat bottom W</b> (flat padded, with seat extension)	<b>Add. charge</b>
Settings for Recaro seat/modular Recaro seating system			
<b>Seat depth SD in cm</b>		<input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46	
(only for Recaro n-Joy-Plus)		<input type="checkbox"/> 48 <input type="checkbox"/> 50 <input type="checkbox"/> 52 <input type="checkbox"/> 54 <input type="checkbox"/> 56	
<b>Lower leg length LL in cm</b>		..... cm (28 ... 54 cm)	
<b>Seat height SH in cm</b>		<input type="checkbox"/> approx. 52 <input type="checkbox"/> approx. 57	
Note: In case of 0° seat inclination, with seat padding, dependent on lower leg length; seat height less lower leg length equals ground clearance, a minimum of 6 cm is recommended; seat height approx. 52 cm not available with mechanical or electric seat tilt			
All Recaro seats with headrest and back bag.			

<b>Seat cushion</b> (only for standard seat)		
EC 06	<input type="checkbox"/>	<b>Seat cushion, black cover</b> (standard equipment) –
EC 10	<input type="checkbox"/>	<b>Seat cushion with incontinence cover</b> Add. charge
EC 07	<input type="checkbox"/>	<b>Without seat cushion</b> Price reduction
EC 11	<input type="checkbox"/>	<b>Terra seat cushion</b> (up to 140 kg/309 lbs) Add. charge
EC 12	<input type="checkbox"/>	<b>Terra Aquos seat cushion</b> (up to 140 kg/309 lbs) Add. charge
EC 13	<input type="checkbox"/>	<b>Aerius seat cushion</b> (up to 140 kg/309 lbs) Add. charge
EC 14	<input type="checkbox"/>	<b>Aerius Active seat cushion</b> (up to 150 kg/330 lbs) Add. charge
EC 17	<input type="checkbox"/>	<b>Terra Flair seat cushion</b> (up to 150 kg/330 lbs) Add. charge
Please refer to our seat cushion brochure for information on additional seat cushions.		

<b>Seat tilt</b>		
EC 50	<input type="checkbox"/>	<b>Electric seat tilt 30° with centre of gravity shifting</b> (seat inclination 0°) Add. charge
EC 54	<input type="checkbox"/>	<b>Electric seat tilt 45° with centre of gravity shifting</b> (seat inclination 0°) Add. charge

<b>Seat height adjustment</b>		
EC 55	<input type="checkbox"/>	<b>Electric seat height adjustment</b> (up to 40 cm, seat inclination 0°) <b>max. load capacity 130 kg/287 lbs!</b> Add. charge

<b>Back adjustment</b>		
ED 20	<input type="checkbox"/>	<b>Mechanical back adjustment 30°</b> (only for standard and contour seat) Add. charge
ED 50	<input type="checkbox"/>	<b>Electric back angle adjustment 30°</b> Add. charge (not available with Recaro n-Joy and Recaro n-Joy-Plus)

<b>Back version</b> (only for standard seats)		
ED 01	<input type="checkbox"/>	<b>Back upholstery</b> (standard equipment) –
ED 02	<input type="checkbox"/>	<b>Without back upholstery</b> Price reduction
Please use our accessories catalogue for the selection of individual back systems!		

<b>Belts</b>		
EP 51	<input type="checkbox"/>	<b>Static four-way chest strap kit</b> (only for Recaro seat/seating system) Add. charge
EP 60	<input type="checkbox"/>	<b>Chest belt with hook-and-loop closure, installed at chest height</b> Add. charge (only for Recaro seat/seating system)
EP 65	<input type="checkbox"/>	<b>Lap belt with buckle</b> (standard equipment) – (Lap belt or four-way chest strap required in Germany for road traffic permit)

<b>Seating system accessories</b>		
EC 20	<input type="checkbox"/>	<b>Abductor, removable</b> (only for standard seat, not available in combination with seat tilt) Add. charge
EC 30	<input type="checkbox"/>	<b>Flip-down abductor</b> (only for standard seat, not available in combination with seat tilt) Add. charge
EE 30	<input type="checkbox"/>	<b>Adapter for headrest installation kit</b> (only for standard and contour seat) Add. charge Please refer to our accessories catalogue for information about headrests.

<b>Footrests</b>		
EB 01	<input type="checkbox"/>	<b>Individual footrests, with plastic footplates</b> (standard equipment) – Available with Standard seat XL and Contour seat XL with max. load capacity of 140 kg/309 lbs; LL up to 50 cm with standard seat. Not available with Standard seat Junior.
EB 02	<input type="checkbox"/>	<b>Individual footrests, with aluminium footplates</b> (continuously angle-adjustable) Add. charge Available with Standard seat XL and Contour seat XL with max. load capacity of 140 kg/309 lbs; LL up to 50 cm with standard seat.
EB 03	<input type="checkbox"/>	<b>Single-panel footrest with aluminium footplate</b> Add. charge (continuously angle-adjustable, LL up to 50 cm with standard seat) Not available with Standard seat Junior.
EB 04	<input type="checkbox"/>	<b>XXL individual footrests with aluminium footplates</b> Add. charge (continuously angle-adjustable, LL up to 50 cm with standard seat) Not available with Standard seat Junior.

<b>EB 07</b>	<input type="checkbox"/>	<b>Without footrest</b>	<b>Price reduction</b>
<b>EB 20</b>	<input type="checkbox"/>	<b>Mechanically elevating footrests with length adjustment and aluminium footplates, angle adjustable</b> (LL with standard seat from 27 cm, with Recaro seat/seating system and contoured seat from 32 cm) Not available with Standard seat XL and Contour seat XL.	<b>Add. charge</b>
<b>EB 50</b>	<input type="checkbox"/>	<b>Electrically elevating footrests with length adjustment and aluminium footplates, jointly controllable</b> (LL from 27 cm for standard seat, from 32 cm for Recaro seat/seating system and contour seat)	<b>Add. charge</b>
<b>EB 51</b>	<input type="checkbox"/>	<b>Individual footrests, electrically elevating, with length adjustment and aluminium footplates</b> (LL from 27 cm for standard seat, from 32 cm for Recaro seat/seating system and contour seat) Not available with Standard seat XL and Contour seat XL.	<b>Add. charge</b>

### Footrest accessories

<b>EB 80</b>	<input type="checkbox"/>	<b>Lateral leg pad for mechanically elevating footrest</b> (pair)	<b>Add. charge</b>
<b>EB 70</b>	<input type="checkbox"/>	<b>Fixation straps for aluminium footplates</b> (included with mounting accessories)	<b>Add. charge</b>

### Armrest

<b>EE 10</b>	<input type="checkbox"/>	<b>Armrest with clothing protector, telescoping</b> (only for standard and contour seat)	-
<b>EE 33</b>	<input type="checkbox"/>	<b>Extended armrest with cut-out</b> (telescoping, only included with Recaro n-joy/njoy+)	-
<b>EE 40</b>	<input type="checkbox"/>	<b>Extended armrest with clothing protector</b> (telescoping, with cut-out for back tube, only for Standard and Contour seat)	<b>Add. charge</b>

### Armrest accessories

<b>EE 34</b>	<input type="checkbox"/>	<b>Channel forearm pad adapter, fixed, mounting side left</b>	<b>Add. charge</b>
<b>EE 24</b>	<input type="checkbox"/>	<b>Channel forearm pad adapter, fixed, mounting side right</b>	<b>Add. charge</b>
<b>EE 35</b>	<input type="checkbox"/>	<b>Channel forearm pad adapter, swivelling, mounting side left</b>	<b>Add. charge</b>
<b>EE 25</b>	<input type="checkbox"/>	<b>Channel forearm pad adapter, swivelling, mounting side right</b>	<b>Add. charge</b>
<b>EE 36</b>	<input type="checkbox"/>	<b>Arm strap, mounting side left</b>	<b>Add. charge</b>
<b>EE 26</b>	<input type="checkbox"/>	<b>Arm strap, mounting side right</b>	<b>Add. charge</b>

For information on channel forearm pads and hand pads, please refer to our accessories catalogue!

### Frame colour

<b>016</b>	<input type="checkbox"/>	<b>Silver metallic</b> (standard equipment)	-
<b>023</b>	<input type="checkbox"/>	<b>Traffic white RAL 9016</b> (standard equipment)	-
<b>024</b>	<input type="checkbox"/>	<b>Traffic blue RAL 5017</b> (standard equipment)	-
<b>022</b>	<input type="checkbox"/>	<b>Cream RAL 9001</b> (standard equipment)	-
<b>038</b>	<input type="checkbox"/>	<b>Sparkle red</b>	<b>Add. charge</b>
<b>021</b>	<input type="checkbox"/>	<b>Sparkle granny-smith</b>	<b>Add. charge</b>
<b>010</b>	<input type="checkbox"/>	<b>Marine blue</b>	<b>Add. charge</b>
<b>020</b>	<input type="checkbox"/>	<b>Saturn blue metallic</b>	<b>Add. charge</b>
<b>005</b>	<input type="checkbox"/>	<b>Sparkle light orange</b>	<b>Add. charge</b>
<b>004</b>	<input type="checkbox"/>	<b>Ice blue</b>	<b>Add. charge</b>
<b>001</b>	<input type="checkbox"/>	<b>Brilliant silver</b>	<b>Add. charge</b>
<b>018</b>	<input type="checkbox"/>	<b>Signal red RAL 3001</b>	<b>Add. charge</b>
<b>002</b>	<input type="checkbox"/>	<b>Anthracite metallic</b>	<b>Add. charge</b>
<b>025</b>	<input type="checkbox"/>	<b>Stainless steel silver</b>	<b>Add. charge</b>
<b>027</b>	<input type="checkbox"/>	<b>Shimano matte</b>	<b>Add. charge</b>
<b>003</b>	<input type="checkbox"/>	<b>Jet black RAL 9005</b>	<b>Add. charge</b>
<b>009</b>	<input type="checkbox"/>	<b>Vinho sparkle</b>	<b>Add. charge</b>
<b>017</b>	<input type="checkbox"/>	<b>Traffic yellow RAL 1023</b>	<b>Add. charge</b>
<b>012</b>	<input type="checkbox"/>	<b>Sky blue</b>	<b>Add. charge</b>
<b>014</b>	<input type="checkbox"/>	<b>Candy red</b>	<b>Add. charge</b>
<b>015</b>	<input type="checkbox"/>	<b>Clear reptile glossy</b>	<b>Add. charge</b>

013	<input type="checkbox"/>	Jet black RAL 9005 matte	Add. charge
019	<input type="checkbox"/>	White pearl	Add. charge
026	<input type="checkbox"/>	Sparkling iron effect light	Add. charge
028	<input type="checkbox"/>	Rainbow	Add. charge

Tread pattern			
EG 01	<input type="checkbox"/>	Ribs	-
EG 02	<input type="checkbox"/>	Lugs	-
EG 50	<input type="checkbox"/>	Caster wheels with ribs, drive wheels with lugs	-

Tyre colour			
EG 03	<input type="checkbox"/>	Grey	-
EG 04	<input type="checkbox"/>	Black (only with lugs)	-

Tyres			
EG 55	<input type="checkbox"/>	Pneumatic tyres (standard equipment)	-
EG 60	<input type="checkbox"/>	Puncture-proof tyres on all wheels	Add. charge
EF 06	<input type="checkbox"/>	Puncture-proof tyres on the caster wheels	Add. charge

Caster wheel diameter			
EF 01	<input type="checkbox"/>	9" caster wheel	-
EF 05	<input type="checkbox"/>	10" caster wheel	-

Caster wheel options			
EF 20	<input type="checkbox"/>	No suspension (standard equipment)	-
EF 02	<input type="checkbox"/>	Suspension	Add. charge

Caster wheel accessories			
EF 03	<input type="checkbox"/>	Mechanical caster swivel lock	Add. charge

Chassis options			
EG 20	<input type="checkbox"/>	No suspension (standard equipment)	-
EG 05	<input type="checkbox"/>	Suspension (drive wheel suspension, incl. drawer for battery case)	Add. charge

Chassis accessories			
ES 05	<input type="checkbox"/>	Electric lighting (required for road traffic permit)	Add. charge
ES 03	<input type="checkbox"/>	Rear marker plate (required for road traffic permit)	Add. charge
EA 10	<input type="checkbox"/>	Splash guard, complete	Add. charge

Drive unit			
ET 01	<input type="checkbox"/>	6 km/h drive unit (standard equipment)	-
ET 09	<input type="checkbox"/>	7.2 km/h/5 mph drive unit (not available for sale in Germany)	-
ET 10	<input type="checkbox"/>	10 km/h drive unit (4-pole-motors) (road traffic permit included, lap belt or four-way chest strap, lighting and rear marker plate required)	Add. charge

Batteries			
EU 03	<input type="checkbox"/>	Without	Price reduction
EU 01	<input type="checkbox"/>	Two 63 Ah (C5) gel batteries (standard equipment)	-
EU 12	<input type="checkbox"/>	Two 63 Ah (C5) AGM batteries	Price reduction

Control unit/control panel			
EU 85	<input type="checkbox"/>	Standard control panel (standard equipment)	
EU 86	<input type="checkbox"/>	Control panel with two external buttons Please indicate button type. <input type="checkbox"/> Piko button (diameter 50 mm, red/black) <input type="checkbox"/> Piko button (diameter 30 mm, red/black)	Add. charge



<b>EU 87</b>	<input type="checkbox"/>	<b>Mid-tray control (MTC) electrically swivelling</b> (charging/programming box required) Please indicate button type. <input type="checkbox"/> <b>Piko button</b> (diameter 50 mm, red/black) <input type="checkbox"/> <b>Piko button</b> (diameter 30 mm, red/black)	<b>Add. charge</b>
--------------	--------------------------	--	--------------------

**Control panel holder**

<b>EU 19</b>	<input type="checkbox"/>	<b>Standard control panel holder</b> (standard equipment)	<b>-</b>
<b>EU 20</b>	<input type="checkbox"/>	<b>Control panel holder, swing-away in parallel</b>	<b>Add. charge</b>
<b>EU 18</b>	<input type="checkbox"/>	<b>Control panel holder, swing-away in parallel, height-adjustable</b>	<b>Add. charge</b>
<b>EU 29</b>	<input type="checkbox"/>	<b>Control panel holder, downwards adjustable</b>	<b>Add. charge</b>

**Control panel holder mounting side** (not for mid-tray control)

	<input type="checkbox"/>	<b>Left</b>	
	<input type="checkbox"/>	<b>Right</b>	

**Tray-side mounting materials for mid-tray control (MTC)**

	<input type="checkbox"/>	<b>Without</b>	
<b>EU 871</b>	<input type="checkbox"/>	<b>Tray with holder</b>	<b>Add. charge</b>
<b>EU 873</b>	<input type="checkbox"/>	<b>Holder only</b>	<b>Add. charge</b>

**External buttons for mid-tray control (MTC)**

	<input type="checkbox"/>	<b>Without</b>	
<b>EU 874</b>	<input type="checkbox"/>	<b>Piko button, diameter 30 mm, red</b>	<b>Add. charge</b>
<b>EU 875</b>	<input type="checkbox"/>	<b>Piko button, diameter 30 mm, yellow</b>	<b>Add. charge</b>
<b>EU 876</b>	<input type="checkbox"/>	<b>Piko button, diameter 30 mm, black</b>	<b>Add. charge</b>
<b>EU 877</b>	<input type="checkbox"/>	<b>Piko button, diameter 50 mm, red</b>	<b>Add. charge</b>
<b>EU 878</b>	<input type="checkbox"/>	<b>Piko button, diameter 50 mm, yellow</b>	<b>Add. charge</b>
<b>EU 879</b>	<input type="checkbox"/>	<b>Piko button, diameter 50 mm, black</b>	<b>Add. charge</b>

One button, "swivel", required; connection of max. 3 buttons, "swivel", "ON/OFF" and "mode" possible.

**Charging/programming box for mid-tray control (MTC)** (includes 1.5 m BUS cable)

<b>EU 56</b>	<input type="checkbox"/>	<b>Yes</b>	<b>Add. charge</b>
--------------	--------------------------	------------	--------------------

**Activation of the drive-away lock**

	<input type="checkbox"/>	<b>Yes</b>	
	<input type="checkbox"/>	<b>No</b>	

**Activation of the electronic steering lock**

	<input type="checkbox"/>	<b>Yes</b>	
	<input type="checkbox"/>	<b>No</b>	

**Control unit accessories**

<b>EU 07</b>	<input type="checkbox"/>	<b>Push-button module to control up to 10 electrical functions</b>	<b>Add. charge</b>
<b>EU 15</b>	<input type="checkbox"/>	<b>Attendant control</b>	<b>Add. charge</b>
<b>EU 22</b>	<input type="checkbox"/>	<b>Protective bow for control panel</b>	<b>Add. charge</b>
<b>EU 23</b>	<input type="checkbox"/>	<b>Fork for tetraplegics, horizontal</b> (adjustable for hand width 70 ... 90 mm)	<b>Add. charge</b>
<b>EU 24</b>	<input type="checkbox"/>	<b>Stick S80</b> (special joystick handle 80 mm)	<b>Add. charge</b>
<b>EU 25</b>	<input type="checkbox"/>	<b>Soft ball</b> (black cover)	<b>Add. charge</b>
<b>EU 26</b>	<input type="checkbox"/>	<b>Ball top, small</b> (diameter 43 mm)	<b>Add. charge</b>
<b>EU 27</b>	<input type="checkbox"/>	<b>Ball top, large</b> (diameter 48 mm)	<b>Add. charge</b>
<b>EU 28</b>	<input type="checkbox"/>	<b>Joystick top, flexible, including ball top, large</b> (Extension by 60 mm; control panel holder, downwards adjustable or swing-away, is recommended)	<b>Add. charge</b>

<b>EU 35</b>	<input type="checkbox"/>	<b>Soft ball, coloured</b> (diameter 30 mm, 4 pieces)	<b>Add. charge</b>
<b>EU 681</b>	<input type="checkbox"/>	<b>2-fold adapter cable with Piko button</b> (diameter 50 mm, black) (free assignment of 2 electrical functions) Please indicate function: <input type="checkbox"/> <b>Footrest, left</b> <input type="checkbox"/> <b>Footrest, right</b> <input type="checkbox"/> <b>Seat tilt</b> <input type="checkbox"/> <b>Back angle adjustment</b> <input type="checkbox"/> <b>Seat height adjustment</b> <input type="checkbox"/> <b>Light</b> <input type="checkbox"/> <b>Warning flasher</b> <input type="checkbox"/> <b>Direction indicator, left</b> <input type="checkbox"/> <b>Direction indicator, right</b> <input type="checkbox"/> <b>Horn</b>	<b>Add. charge</b>
<b>EU 682</b>	<input type="checkbox"/>	<b>2-fold adapter cable with Piko button</b> (diameter 30 mm, black) (free assignment of 2 electrical functions) Please indicate function: <input type="checkbox"/> <b>Footrest, left</b> <input type="checkbox"/> <b>Footrest, right</b> <input type="checkbox"/> <b>Seat tilt</b> <input type="checkbox"/> <b>Back angle adjustment</b> <input type="checkbox"/> <b>Seat height adjustment</b> <input type="checkbox"/> <b>Light</b> <input type="checkbox"/> <b>Warning flasher</b> <input type="checkbox"/> <b>Direction indicator, left</b> <input type="checkbox"/> <b>Direction indicator, right</b> <input type="checkbox"/> <b>Horn</b>	<b>Add. charge</b>

<b>General accessories</b>			
<b>ES 20</b>	<input type="checkbox"/>	<b>Pannier bag, detachable</b>	<b>Add. charge</b>
<b>ES 22</b>	<input type="checkbox"/>	<b>Accessory pouch</b> (e.g. for battery charger, tool kit)	<b>Add. charge</b>
<b>ES 24</b>	<input type="checkbox"/>	<b>Pocket for mobile phone</b>	<b>Add. charge</b>
<b>ES 32</b>	<input type="checkbox"/>	<b>Cover cap for standard seating units</b>	<b>Add. charge</b>
<b>ES 35</b>	<input type="checkbox"/>	<b>Tray, swing-away to the side</b> (includes installation accessories; cut-out for control panel must be adapted)	<b>Add. charge</b>
<b>EZ 80</b>	<input type="checkbox"/>	<b>Airman pump</b>	<b>Add. charge</b>
<b>EA 62</b>	<input type="checkbox"/>	<b>Rear bumper</b>	<b>Add. charge</b>
<b>ES 37</b>	<input type="checkbox"/>	<b>External odometer</b> (only in combination with protective bow for control panel)	<b>Add. charge</b>
<b>EA 60</b>	<input type="checkbox"/>	<b>Crutch holder</b>	<b>Add. charge</b>
<b>ES 38</b>	<input type="checkbox"/>	<b>Folding rearview mirror, mounting side left</b>	<b>Add. charge</b>
<b>ES 28</b>	<input type="checkbox"/>	<b>Folding rearview mirror, mounting side right</b>	<b>Add. charge</b>
(Mounting to the control panel side not possible with holder, swing-away in parallel; not available in combination with holder, swing-away in parallel, and tray)			
<b>EB 52</b>	<input type="checkbox"/>	<b>Tool kit</b>	<b>Add. charge</b>
<b>ES 60</b>	<input type="checkbox"/>	<b>Vehicle transportation kit</b> (required for AMF anchor point system) (only possible for user weight of 100 kg max./220 lbs) (only possible with standard and contour seat)	<b>Add. charge</b>

<b>Special options</b>			
<b>24817.001</b>	<input type="checkbox"/>	<b>Special colour two-layer coating</b>	<b>Add. charge</b>
<b>24817.000</b>	<input type="checkbox"/>	<b>Special colour single-layer coating</b>	<b>Add. charge</b>
<b>21187.002</b>	<input type="checkbox"/>	<b>Special dimensions for platform for seating shell</b>	<b>Add. charge</b>
<b>20465.016</b>	<input type="checkbox"/>	<b>Standard back tube extension/shortening</b>	<b>Add. charge</b>
<b>40405.002</b>	<input type="checkbox"/>	<b>Back upholstery for power wheelchairs with special length</b>	<b>Add. charge</b>
<b>20465.011</b>	<input type="checkbox"/>	<b>Contour back with standard seat bottom or contour seat bottom with standard back</b>	<b>Add. charge</b>
<b>36362.001</b>	<input type="checkbox"/>	<b>Avantgarde2 footplates, locked</b>	<b>Add. charge</b>
<b>36362.000</b>	<input type="checkbox"/>	<b>Avantgarde2 footplates, unlocked</b>	<b>Add. charge</b>
<b>25142.001</b>	<input type="checkbox"/>	<b>Amputation legrest for power wheelchairs</b>	<b>Add. charge</b>



<b>22298.006</b>	<input type="checkbox"/>	<b>Removable enAble control panel for attendant</b>	<b>Add. charge</b>
<b>22517.005</b>	<input type="checkbox"/>	<b>"ALU" e-Hockey sport bumper</b> (available only in combination with special option 22517.003 and 59060.000)	<b>Add. charge</b>
<b>22517.003</b>	<input type="checkbox"/>	<b>e-Hockey ball protection side</b> (available only in combination with special option 22517.005 and 59060.000)	<b>Add. charge</b>
<b>59060.000</b>	<input type="checkbox"/>	<b>e-Hockey hockey stick</b> (available only in combination with special option 22517.005 and 22517.003)	<b>Add. charge</b>