

B600 Power wheelchair

Max. load capacity: 140 kg/309 lbs
(standard equipment)

Article number: 490E32=10000_K



Standard Equipment

- 6 km/h
- enAble50 control unit
- 63 Ah (C5) gel batteries
- Battery charger
- Lights
- Standard seat
- Seat cushion
- Drive wheel suspension
- Armrests with clothing protector
- Footrests
- Bumper
- Lap belt with buckle

Mandatory information*

Body height in cm _____

Body weight in kg _____

1 Seat depth _____

2 Trunk depth _____

3 Armrest height _____

4 Armpit height _____

5 Shoulder height _____

6 Head height _____

7 Forearm length to the middle of the palm _____

8 Trunk width _____

9 Lower leg length _____

10 Seat width _____


Order no.: _____	
-------------------------	--

<input type="checkbox"/> Provide estimate only <input type="checkbox"/> Order	
<input type="checkbox"/> Attachment for custom fabrication enclosed	
Billing address/customer number	
Delivery address	
Sales person	Name Date Signature

Basic data (voluntary information)

Gender _____ Indication _____

Age _____ Paying party _____



* This information will help to support user-oriented settings and product development. Ottobock guarantees the confidential and purely internal use of these data.

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please use the spare parts catalogue for spare parts orders.

Information



EC 01/02 Standard seat



EC 03/04 Contour seat



EC 40 Recaro N-joy



EC 42 Backrest LX



EC 43 Backrest LT



EU 07 Push-button module



EC 44 Seat bottom F



EC 45 Seat bottom X



EC 46 Seat bottom W



EC 50 Electric seat tilt



Electric back angle
adjustment 30°



EC 55 Lift seat



EB 01 Individual footrests



EB 03 Single-panel footrest



EB 20 Elevating footrest

Standard seat			
EC 18	<input type="checkbox"/>	Standard seat, Junior (SD 34 ... 40 cm, SW 34 ... 40 cm)	Add. charge
EC 01	<input type="checkbox"/>	Standard seat, small (SD 38 ... 46 cm, SW 38 ... 42 cm)	-
EC 02	<input type="checkbox"/>	Standard seat, large (SD 42 ... 50 cm, SW 43 ... 48 cm)	-
EC 19	<input type="checkbox"/>	Standard seat, XL (SD 42 ... 50 cm, SW 50 ... 56 cm)	Add. charge
Note: Standard seat, XL is not in combination with SH 45 cm possible. Standard seat, XL is not in combination with footrests possible: Optioncode: EB 01, EB 02, EB 20, EB 51.			
Settings for the standard seat			
Back angle		<input type="checkbox"/> -9/1/11/21° <input type="checkbox"/> 0/10/20/30°	
Note: Not available with mechanical and electric back adjustment. Back angle -9/1/11/21° is not possible in combination with SD 34 and 36 cm.			
Seat inclination SI		<input type="checkbox"/> 0° <input type="checkbox"/> 3° <input type="checkbox"/> 6°	
EC 65	<input type="checkbox"/>	Seat inclination 20° (only for standard seat, SH 45 and 50 cm)	Add. charge
3°, 6°, and 20° seat inclination not available with electric seat tilt and seat height adjustment.			
20° seat inclination not available with 180 kg/397 lbs chassis!			
Seat depth SD in cm		<input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50	
Seat width SW	 cm	
Lower leg length LL	 cm (25 ... 54 cm) without seat cushion	
Seat height SH in cm		<input type="checkbox"/> 45 <input type="checkbox"/> 50 <input type="checkbox"/> 55	
(Note: for 3° seat inclination, without seat cushion, dependent on LL, not available with electric seat height adjustment) (SH minus LL equals ground clearance, minimum of 6 cm is recommended)			
Contour seat			
EC 03	<input type="checkbox"/>	Contour seat, small (SD 36 ... 44 cm, SW 38 ... 42 cm)	Add. charge
EC 04	<input type="checkbox"/>	Contour seat, large (SD 40 ... 48 cm, SW 43 ... 48 cm)	Add. charge
EC 21	<input type="checkbox"/>	Contour seat, XL (SD 40 ... 48 cm, SW 50 ... 56 cm) (not possible in combination with seat height 52 cm)	Add. charge
Settings for the contour seat			
Seat depth SD in cm		<input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48	
Seat width SW	 cm	
Lower leg length LL	 cm (28 ... 54 cm)	
Seat height SH in cm		<input type="checkbox"/> approx. 52 <input type="checkbox"/> approx. 57	
(Note: for seat inclination of 0°, with seat padding, dependent on LL, not available with electric seat height adjustment) (SH minus LL equals ground clearance, minimum of 6 cm is recommended)			
Recaro seat			
EC 40	<input type="checkbox"/>	Recaro n-Joy (manually adjustable backrest)	Add. charge
EC 41	<input type="checkbox"/>	Recaro n-Joy-Plus (adaptable on site, manually adjustable backrest) Seat depth SD 38 ... 56 cm	Add. charge
Modular Recaro seating system			
EC 42	<input type="checkbox"/>	Backrest LX (manually adjustable)	Add. charge
EC 43	<input type="checkbox"/>	Backrest LT (lateral supports with different setting possibilities, manually adjustable)	Add. charge
EC 44	<input type="checkbox"/>	Seat bottom F (flat padded, seat bottom sewn to lateral guides)	Add. charge
EC 45	<input type="checkbox"/>	Seat bottom X (contoured lateral guides with seat extension)	Add. charge
EC 46	<input type="checkbox"/>	Seat bottom W (flat padded, with seat extension)	Add. charge
Settings for Recaro seat/modular Recaro seating system			
Seat depth SD in cm		<input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46	
(only for Recaro-n-Joy-Plus)		<input type="checkbox"/> 48 <input type="checkbox"/> 50 <input type="checkbox"/> 52 <input type="checkbox"/> 54 <input type="checkbox"/> 56	
Lower leg length LL	 cm (28 ... 54 cm)	
Seat height SH		<input type="checkbox"/> approx. 52 cm <input type="checkbox"/> approx. 57 cm	
(Note: for seat inclination of 0°, with seat padding, not available with electric seat height adjustment), SH minus LL equals ground clearance, minimum of 6 cm is recommended, SH approx. 52 cm not available with electric seat tilt and 180 kg/397 lbs chassis!)			
All Recaro seats include headrest and back bag.			

Chassis options		
EA 14	<input type="checkbox"/>	140 kg/309 lbs chassis (standard equipment) -
EA 15	<input type="checkbox"/>	180 kg/397 lbs chassis incl. single-panel aluminium footplate and puncture-proof tyres (enhanced drive wheel suspension, reinforced side panels and seat frame, 10 km/h version incl. more powerful 10 km/h drive unit, wheelbase 48 ... 56 cm) Add. charge

Telescoping wheelbase		
	<input type="checkbox"/> cm 44 ... 56 cm, for 180 kg/397 lbs chassis 48 ... 56 cm - (depending on caster wheel options, standard setting 46 cm)

Chassis accessories		
ES 05	<input type="checkbox"/>	Rear marker plate (required in Germany for road traffic permit) Add. charge
EA 10	<input type="checkbox"/>	Splash guard, complete Add. charge

Electric lighting		
ES 07	<input type="checkbox"/>	No lighting Price reduction

Seat cushion (only for standard seat)		
EC 06	<input type="checkbox"/>	Seat cushion, black cover (standard equipment) -
EC 10	<input type="checkbox"/>	Seat cushion with incontinence cover Add. charge
EC 07	<input type="checkbox"/>	Without seat cushion Price reduction
EC 11	<input type="checkbox"/>	Terra seat cushion (up to 140 kg/309 lbs) Add. charge
EC 12	<input type="checkbox"/>	Terra Aquos seat cushion (up to 140 kg/309 lbs) Add. charge
EC 13	<input type="checkbox"/>	Aerius seat cushion (up to 140 kg/309 lbs) Add. charge
EC 14	<input type="checkbox"/>	Aerius Active seat cushion (up to 150 kg/330 lbs) Add. charge
EC 17	<input type="checkbox"/>	Terra Flair seat cushion (up to 150 kg/330 lbs) Add. charge

Please refer to our seat cushion brochure for information on additional seat cushions.

Seat tilt		
EC 50	<input type="checkbox"/>	Electric seat tilt 30° with centre of gravity shifting (seat inclination 0°) Add. charge
EC 54	<input type="checkbox"/>	Electric seat tilt 45° with centre of gravity shifting (seat inclination 0°) Add. charge

Seat height adjustment		
EC 55	<input type="checkbox"/>	Electric seat height adjustment (up to 40 cm, seat inclination 0°) Add. charge (Available only in combination with 180 kg chassis! Max. load capacity: 180 kg/397 lbs!)
EC 60	<input type="checkbox"/>	Electric seat height adjustment "Heavy Duty" (up to 30 cm, seat inclination 0°) Add. charge max. load capacity 180 kg/397 lbs!

Back adjustment		
ED 50	<input type="checkbox"/>	Electric back angle adjustment 30° Add. charge (not available with Recaro n-Joy and Recaro n-Joy-Plus)

Back version (only for standard seat)		
ED 01	<input type="checkbox"/>	Back upholstery (standard equipment) -
ED 02	<input type="checkbox"/>	Without back upholstery -

Belts		
EP 51	<input type="checkbox"/>	Static four-way chest strap kit (only for Recaro seat/seating system) Add. charge (not available in combination with electric seat tilt and electric seat height adjustment)
EP 60	<input type="checkbox"/>	Chest belt with hook-and-loop closure, installed at chest height (only for Recaro seat/seating system) Add. charge
EP 65	<input type="checkbox"/>	Lap belt with buckle (standard equipment) Add. charge

Lap belt or four-way chest strap required in Germany for road traffic permit.

Seating system accessories			
EC 20	<input type="checkbox"/>	Abductor, detachable (only for standard seat, cannot be combined with electric seat tilt and seat height adjustment)	Add. charge
EC 30	<input type="checkbox"/>	Flip-down abductor (only for standard seat, not available in combination with electric seat tilt and seat height adjustment)	Add. charge
EE 30	<input type="checkbox"/>	Adapter for headrest installation kit (only for Standard and Contour seat)	Add. charge

Please refer to our accessories catalogue for information about headrests.

Footrests			
EB 01	<input type="checkbox"/>	Individual footrests, plastic footplates (standard equipment) Available with Standard seat XL and Contour seat XL with max. load capacity of 140 kg/309 lbs; LL up to 50 cm with standard seat. Not available with Standard seat Junior.	–
EB 02	<input type="checkbox"/>	Individual footrests, with aluminium footplates (continuously angle-adjustable) Available with Standard seat XL and Contour seat XL with max. load capacity of 140 kg/309 lbs; LL up to 50 cm with standard seat.	Add. charge
EB 03	<input type="checkbox"/>	Single-panel footrest with aluminium footplate (continuously angle-adjustable, LL up to 50 cm with standard seat) Not available with Standard seat Junior.	Add. charge
EB 04	<input type="checkbox"/>	XXL individual footrests with aluminium footplates (continuously angle-adjustable, LL up to 50 cm with standard seat) Not available with Standard seat Junior.	Add. charge
EB 07	<input type="checkbox"/>	Without footrest	Price reduction
EB 20	<input type="checkbox"/>	Mechanically elevating footrests with length adjustment and aluminium footplates angle-adjustable, max. load capacity 140 kg/309 lbs! (LL from 27 cm for Standard seat, from 32 cm for Recaro seat/seating system and Contour seat) Not available with Standard seat XL and Contour seat XL.	Add. charge
EB 50	<input type="checkbox"/>	Electrically elevating footrest with length adjustment and aluminium footplates (up to 140 kg/309 lbs incl. individual aluminium footplates, over 140 kg incl. single-panel aluminium footplate; LL from 27 cm for Standard seat, from 32 cm for Recaro seat/seating system and Contour seat) Not available with Standard seat Junior.	Add. charge

Footrest accessories			
EB 80	<input type="checkbox"/>	Lateral leg pad for electrically and mechanically elevating footrest (pair)	Add. charge
EB 70	<input type="checkbox"/>	Fixation straps for aluminium footplates (included with installation accessories)	Add. charge

Armrest			
EE 10	<input type="checkbox"/>	Armrest with clothing protector (standard equipment, only for Standard and Contour seats)	–
EE 33	<input type="checkbox"/>	Extended armrest with cut-out (telescoping, only included with Recaro n-joy/njoy+)	–
EE 40	<input type="checkbox"/>	Extended armrest with clothing protector (telescoping, with cut-out for back tube, only for Standard and Contour seat)	Add. charge

Armrest accessories			
EE 34	<input type="checkbox"/>	Channel forearm pad adapter, fixed, mounting side left	Add. charge
EE 24	<input type="checkbox"/>	Channel forearm pad adapter, fixed, mounting side right	Add. charge
EE 35	<input type="checkbox"/>	Channel forearm pad adapter, swivelling, mounting side left	Add. charge
EE 25	<input type="checkbox"/>	Channel forearm pad adapter, swivelling, mounting side right	Add. charge
EE 36	<input type="checkbox"/>	Arm strap, mounting side left	Add. charge
EE 26	<input type="checkbox"/>	Arm strap, mounting side right	Add. charge

For information on channel forearm pads and hand pads, please refer to our accessories catalogue!

Frame colour			
016	<input type="checkbox"/>	Silver metallic (standard equipment)	–
023	<input type="checkbox"/>	Traffic white RAL 9016 (standard equipment)	–
024	<input type="checkbox"/>	Traffic blue RAL 5017 (standard equipment)	–
022	<input type="checkbox"/>	Cream RAL 9001 (standard equipment)	–
038	<input type="checkbox"/>	Sparkle red	Add. charge
021	<input type="checkbox"/>	Sparkle granny-smith	Add. charge

010	<input type="checkbox"/>	Marine blue	Add. charge
020	<input type="checkbox"/>	Saturn blue metallic	Add. charge
005	<input type="checkbox"/>	Sparkle light orange	Add. charge
004	<input type="checkbox"/>	Ice blue	Add. charge
001	<input type="checkbox"/>	Brilliant silver	Add. charge
018	<input type="checkbox"/>	Signal red RAL 3001	Add. charge
002	<input type="checkbox"/>	Anthracite metallic	Add. charge
025	<input type="checkbox"/>	Stainless steel silver	Add. charge
027	<input type="checkbox"/>	Shimano matte	Add. charge
003	<input type="checkbox"/>	Jet black RAL 9005	Add. charge
009	<input type="checkbox"/>	Vinho sparkle	Add. charge
017	<input type="checkbox"/>	Traffic yellow RAL 1023	Add. charge
012	<input type="checkbox"/>	Sky blue	Add. charge
014	<input type="checkbox"/>	Candy red	Add. charge
015	<input type="checkbox"/>	Clear reptile glossy	Add. charge
013	<input type="checkbox"/>	Jet black RAL 9005 matte	Add. charge
019	<input type="checkbox"/>	White pearl	Add. charge
026	<input type="checkbox"/>	Sparkling Iron Effect Light	Add. charge
028	<input type="checkbox"/>	Rainbow	Add. charge

Tread pattern			
EG 01	<input type="checkbox"/>	Ribs	-
EG 02	<input type="checkbox"/>	Lugs	-
EG 50	<input type="checkbox"/>	Caster wheels with ribs, drive wheels with lugs	-

Tyre colour			
EG 03	<input type="checkbox"/>	Grey	-
EG 04	<input type="checkbox"/>	Black (only with lugs)	-

Tyres			
EG 55	<input type="checkbox"/>	Pneumatic tyres (standard equipment)	-
EG 60	<input type="checkbox"/>	Puncture-proof tyres on all wheels	Add. charge
EF 06	<input type="checkbox"/>	Puncture-proof tyres on the caster wheels	Add. charge

Caster wheel diameter			
EF 01	<input type="checkbox"/>	9" caster wheel	-
EF 05	<input type="checkbox"/>	10" caster wheel	-

Caster wheel options			
EF 20	<input type="checkbox"/>	No suspension (standard equipment)	-
EF 02	<input type="checkbox"/>	Suspension	Add. charge

Caster wheel accessories			
EF 03	<input type="checkbox"/>	Mechanical caster swivel lock	Add. charge
EF 04	<input type="checkbox"/>	Curb climbing assist (not available in combination with single-panel footplate and electric seat height adjustment, only in combination with 46 and 48 cm wheelbase)	Add. charge

Drive unit			
ET 01	<input type="checkbox"/>	6 km/h drive unit (standard equipment)	-
ET 10	<input type="checkbox"/>	10 km/h drive unit (4-pole-motors, road traffic permit included, lap belt or four-way chest strap, lighting and rear marker plate required)	Add. charge

Batteries		
EU 03	<input type="checkbox"/>	Without Price reduction
EU 01	<input type="checkbox"/>	Two 63 Ah (C5) gel batteries (standard equipment) -
EU 12	<input type="checkbox"/>	Two 63 Ah (C5) AGM batteries Price reduction

Battery charger		
EU 91	<input type="checkbox"/>	Battery charger (standard equipment) -
EU 02	<input type="checkbox"/>	Battery charger without fan (low-noise) for gel batteries Add. charge

Control unit/control panel		
EU 85	<input type="checkbox"/>	Standard control panel (standard equipment) -
EU 86	<input type="checkbox"/>	Control panel with two external buttons Please indicate button type! <input type="checkbox"/> Piko button (diameter 50 mm, red/black) <input type="checkbox"/> Piko button (diameter 30 mm, red/black) Add. charge
EU 87	<input type="checkbox"/>	Mid-tray control (MTC) electrically swivelling (charging/programming box required) Please indicate button type! <input type="checkbox"/> Piko button (diameter 50 mm, red/black) <input type="checkbox"/> Piko button (diameter 30 mm, red/black) Add. charge

Control panel holder		
EU 19	<input type="checkbox"/>	Standard control panel holder (standard equipment) -
EU 20	<input type="checkbox"/>	Control panel holder, swing-away in parallel Add. charge
EU 18	<input type="checkbox"/>	Control panel holder, swing-away, height-adjustable Add. charge
EU 29	<input type="checkbox"/>	Control panel holder, downwards adjustable Add. charge

Control panel holder mounting side (not for mid-tray control)		
	<input type="checkbox"/>	Right
	<input type="checkbox"/>	Left

Tray-side mounting materials for mid-tray control (MTC)		
	<input type="checkbox"/>	Without
EU 871	<input type="checkbox"/>	Tray with holder Add. charge
EU 873	<input type="checkbox"/>	Holder only Add. charge

External buttons for mid-tray control (MTC)		
	<input type="checkbox"/>	Without
EU 874	<input type="checkbox"/>	Piko button, diameter 30 mm, red Add. charge
EU 875	<input type="checkbox"/>	Piko button, diameter 30 mm, yellow Add. charge
EU 876	<input type="checkbox"/>	Piko button, diameter 30 mm, black Add. charge
EU 877	<input type="checkbox"/>	Piko button, diameter 50 mm, red Add. charge
EU 878	<input type="checkbox"/>	Piko button, diameter 50 mm, yellow Add. charge
EU 879	<input type="checkbox"/>	Piko button, diameter 50 mm, black Add. charge

One button, "swivel", required; connection of max. 3 buttons, "swivel", "ON/OFF" and "mode" possible.

Charging/programming box for mid-tray control (MTC) (includes 1.5 m BUS cable)		
EU 56	<input type="checkbox"/>	Yes Add. charge

Special controls		
EU 07	<input type="checkbox"/>	Push-button module (control of max. 10 electric options without LCD) Add. charge
EU 97	<input type="checkbox"/>	Programmable LCD, IR environment control, Bluetooth and LED control included Add. charge
EU 98	<input type="checkbox"/>	Programmable LCD, IR environment control and LED control included Add. charge

The following special controls require an LCD:		
EU 09	<input type="checkbox"/>	Chin control with mini joystick, electrically swing-away (standard control panel included) Please indicate: <input type="checkbox"/> Soft ball with elastomer <input type="checkbox"/> Soft ball with spring, hard <input type="checkbox"/> Soft ball with spring, soft Please indicate installation version: <input type="checkbox"/> Chin swivel arm, left <input type="checkbox"/> Chin swivel arm, right <input type="checkbox"/> Chin swivel arm overhead, mounting side left (multifunction switch MFS right) <input type="checkbox"/> Chin swivel arm overhead, mounting side right (multifunction switch MFS left) Add. charge
	<input type="checkbox"/>	Chin control with midi joystick, electrically swing-away (standard control panel included) Please indicate installation version: <input type="checkbox"/> Chin swivel arm, left <input type="checkbox"/> Chin swivel arm, right <input type="checkbox"/> Chin swivel arm overhead, mounting side left (multifunction switch MFS right) <input type="checkbox"/> Chin swivel arm overhead, mounting side right (multifunction switch MFS left) Add. charge
EU 10	<input type="checkbox"/>	Foot control, left (standard control panel included) Add. charge
EU 10	<input type="checkbox"/>	Foot control, right (standard control panel included) Add. charge
EU 11	<input type="checkbox"/>	Sip and puff control (only with 6 km/h drive unit and Recaro seat/seating system; standard control panel included) Add. charge
EU 40	<input type="checkbox"/>	Calibratable joystick, standard, includes multifunction button (standard control panel included) Add. charge
EU 45	<input type="checkbox"/>	Calibratable joystick, includes multifunction button, for assembly in tray (includes holder for LCD with gooseneck for attachment in the tray and standard control panel) Add. charge
Joystick top version: Please indicate.		
	<input type="checkbox"/>	Control knob with elastomer
	<input type="checkbox"/>	Soft ball with elastomer
	<input type="checkbox"/>	Soft ball with spring, hard
	<input type="checkbox"/>	Soft ball with spring, soft
EU 68	<input type="checkbox"/>	Adapter for 1 – 6 buttons (required for 1 – 6 button control; select buttons under U 691 to 704) Add. charge
The following options require the IR LCD display code EU 97 or EU 98:		
EU 61	<input type="checkbox"/>	IR remote control hands-free telephone including mains adapter Add. charge
EU 63	<input type="checkbox"/>	“Dragger 32” mouse click input assistance Add. charge
EU 59	<input type="checkbox"/>	Adapter cable for left/right mouse click Add. charge <input type="checkbox"/> Piko button (diameter 50 mm) (red/black) <input type="checkbox"/> Piko button (diameter 30 mm) (red/black)
Special control accessories		
EU 69	<input type="checkbox"/>	Easywave wireless module Add. charge
EU 67	<input type="checkbox"/>	ON/OFF or EMERGENCY STOP with gooseneck for LCD Add. charge (ON/OFF button for LED control required; select button under EU 691 to EU 704)
EU 75	<input type="checkbox"/>	ON/OFF button for LCD module, includes MFS adapter cable Add. charge (ON/OFF button for LED control required; select button under EU 691 to EU 704)
Piko button (diameter 50 mm)		
EU 691	<input type="checkbox"/>	Red Add. charge
EU 693	<input type="checkbox"/>	Black Add. charge
EU 692	<input type="checkbox"/>	Yellow Add. charge
EU 697	<input type="checkbox"/>	Silver Add. charge

Piko button (diameter 30 mm)			
EU 701	<input type="checkbox"/>	Red	Add. charge
EU 702	<input type="checkbox"/>	Yellow	Add. charge
EU 703	<input type="checkbox"/>	Black	Add. charge
EU 704	<input type="checkbox"/>	Silver	Add. charge
Activation of the drive-away lock			
	<input type="checkbox"/>	Yes	
	<input type="checkbox"/>	No	
Activation of the electronic steering lock			
	<input type="checkbox"/>	Yes	
	<input type="checkbox"/>	No	
LCD holder			
EU 31	<input type="checkbox"/>	Standard LCD holder (standard equipment for LCD)	-
EU 30	<input type="checkbox"/>	LCD holder, swing-away, height- and angle-adjustable	Add. charge
EU 32	<input type="checkbox"/>	LCD holder with gooseneck	Add. charge
EU 34	<input type="checkbox"/>	LCD holder with gooseneck for tray (included)	Add. charge
Control unit accessories			
EU 15	<input type="checkbox"/>	Attendant control	Add. charge
EU 22	<input type="checkbox"/>	Protective bow for standard control panel	Add. charge
EU 23	<input type="checkbox"/>	Fork for tetraplegics, horizontal (adjustable for hand width 70 ... 90 mm)	Add. charge
EU 24	<input type="checkbox"/>	Stick S80 (special joystick handle 80 mm)	Add. charge
EU 25	<input type="checkbox"/>	Soft ball (black cover)	Add. charge
EU 26	<input type="checkbox"/>	Ball top, small (diameter 43 mm)	Add. charge
EU 27	<input type="checkbox"/>	Ball top, large (diameter 48 mm)	Add. charge
EU 28	<input type="checkbox"/>	Joystick top, flexible, including ball top, large	Add. charge
Extension by 60 mm; control panel holder, downwards adjustable or swing-away, is recommended.			
EU 35	<input type="checkbox"/>	Soft ball, coloured (diameter 30 mm, 4 pieces)	Add. charge
EU 50	<input type="checkbox"/>	Handheld programming device (dealer version, not subject to discount)	Add. charge
EU 55	<input type="checkbox"/>	Programming software enAble50 with USB adapter (not subject to discount, certificate required!)	Add. charge
EU 681	<input type="checkbox"/>	2-fold adapter cable with Piko button (diameter 50 mm, black) (free assignment of 2 electrical functions) Please indicate function: <input type="checkbox"/> Footrest, left <input type="checkbox"/> Footrest, right <input type="checkbox"/> Seat tilt <input type="checkbox"/> Back angle adjustment <input type="checkbox"/> Seat height adjustment <input type="checkbox"/> Light <input type="checkbox"/> Warning flasher <input type="checkbox"/> Direction indicator, left <input type="checkbox"/> Direction indicator, right <input type="checkbox"/> Horn	Add. charge

EU 682	<input type="checkbox"/>	2-fold adapter cable with Piko button (diameter 30 mm, black) (free assignment of 2 electrical functions) Please indicate function: <input type="checkbox"/> Footrest, left <input type="checkbox"/> Footrest, right <input type="checkbox"/> Seat tilt <input type="checkbox"/> Back angle adjustment <input type="checkbox"/> Seat height adjustment <input type="checkbox"/> Light <input type="checkbox"/> Warning flasher <input type="checkbox"/> Direction indicator, left <input type="checkbox"/> Direction indicator, right <input type="checkbox"/> Horn	Add. charge
---------------	--------------------------	--	--------------------

General accessories			
ES 20	<input type="checkbox"/>	Pannier bag, detachable	Add. charge
ES 22	<input type="checkbox"/>	Accessory pouch (e.g. for battery charger, tool kit)	Add. charge
ES 24	<input type="checkbox"/>	Pocket for mobile phone	Add. charge
ES 32	<input type="checkbox"/>	Cover cap for standard seating units	Add. charge
ES 35	<input type="checkbox"/>	Tray, swing-away to the side (incl. mounting accessories; cut-out for control panel must be adapted)	Add. charge
EZ 80	<input type="checkbox"/>	Airman pump	Add. charge
EA 60	<input type="checkbox"/>	Crutch holder	Add. charge
ES 38	<input type="checkbox"/>	Folding rearview mirror, mounting side left	Add. charge
ES 28	<input type="checkbox"/>	Folding rearview mirror, mounting side right (Mounting to the control panel side not possible with holder, swing-away in parallel; not available in combination with holder, swing-away in parallel, and tray)	Add. charge
ES 37	<input type="checkbox"/>	External kilometre indicator (only in combination with protective bow for control panel)	Add. charge
EZ 52	<input type="checkbox"/>	Tool kit	Add. charge
ES 60	<input type="checkbox"/>	Vehicle transportation kit (required for AMF anchor point system, only available in combination with max. load capacity 100 kg/220 lbs)	Add. charge
ES 10	<input type="checkbox"/>	Platform for respirator	Add. charge

Special options			
24817.001	<input type="checkbox"/>	Special colour two-layer coating	Add. charge
24817.000	<input type="checkbox"/>	Special colour single-layer coating	Add. charge
21187.002	<input type="checkbox"/>	Special dimensions for platform for seating shell	Add. charge
20465.016	<input type="checkbox"/>	Standard back tube extension/shortening	Add. charge
40405.002	<input type="checkbox"/>	Back upholstery for power wheelchairs with special length	Add. charge
20465.011	<input type="checkbox"/>	Contour back with standard seat bottom or contour seat bottom with standard back	Add. charge
36362.001	<input type="checkbox"/>	Avantgarde2 footplates, locked	Add. charge
36362.000	<input type="checkbox"/>	Avantgarde2 footplates, unlocked	Add. charge
22298.006	<input type="checkbox"/>	Removable enAble control panel for attendant	Add. charge
22298.007	<input type="checkbox"/>	Removable enAble control panel for attendant, angle-adjustable	Add. charge
49058.001	<input type="checkbox"/>	Stand-alone joystick, replaces control panel	Add. charge
26156.000	<input type="checkbox"/>	Replace satellite switch, by 2 buttons (chin control)	Add. charge
22569.004	<input type="checkbox"/>	Lip joystick (chin control)	Add. charge
22762.003	<input type="checkbox"/>	Push-button module and control panel on the same side	Add. charge
19893.001	<input type="checkbox"/>	E-series wireless stop	Add. charge
23279.000	<input type="checkbox"/>	Mid-tray control according to drawing	Add. charge
23279.001	<input type="checkbox"/>	Mid-tray control for electrically swing-away control panel	Add. charge
22517.005	<input type="checkbox"/>	"ALU" e-Hockey sport bumper (available only in combination with special option 22517.003 and 59060.000)	Add. charge

22517.003	<input type="checkbox"/>	e-Hockey ball protection side (available only in combination with special option 22517.005 and 59060.000)	Add. charge
59060.000	<input type="checkbox"/>	e-Hockey hockey stick (available only in combination with special option 22517.005 and 22517.003)	Add. charge
25142.001	<input type="checkbox"/>	Amputation legrest for power wheelchairs	Add. charge