

Order no.:

C1000 DS

Power wheelchair **E101**

Max. load capacity: 140 kg/309 lbs
(standard equipment)

Article number: 490E68=00000_K

C1000 DS/S

Power wheelchair **E102**

Max. load capacity: 140 kg/309 lbs
(standard equipment)

C1000 DS/H

Power wheelchair **E103**

Max. load capacity: 130 kg/287 lbs
(standard equipment)

Provide estimate only

Attachment for custom fabrication enclosed

Billing address/customer number

Delivery address

Sales person

Name

.....
Date

Signature



Figure shows add. charge options

Standard Equipment


- enAble 50 impulse control
- 6 km/h
- 63 Ah (C5) gel batteries
- Battery charger
- Lighting
- Luggage carrier
- Suspension on all wheels
- Standard seat
- Seat cushion
- Armrests
- Individual footrests
(up to 140 kg/308 lbs load capacity)
- Splash guard, complete
- Lap belt with buckle
- Chassis colour Sparkling Iron Effect Light
- With C1000 DS/H, incl. electric seat height adjustment (max. load capacity 130 kg/286 lbs)

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please use the spare parts catalogue for spare parts orders.

Basic data (voluntary information)

Gender _____ Indication _____

Age _____ Paying party _____



Mandatory information*

Body height in cm _____

Body weight in kg _____

1 Seat depth _____

2 Trunk depth _____

3 Armrest height _____

4 Armpit height _____

5 Shoulder height _____

6 Head height _____

7 Forearm length to the middle of the palm _____

8 Trunk width _____

9 Lower leg length _____

10 Seat width _____

* This information will help to support user-oriented settings and product development. Ottobock guarantees the confidential and purely internal use of these data.

Standard seat		
EC 01	<input type="checkbox"/>	Standard seat, small (SD 38 ... 46 cm, SW 38 ... 42 cm) -
EC 02	<input type="checkbox"/>	Standard seat, large (SD 42 ... 50 cm, SW 43 ... 48 cm) -
EC 19	<input type="checkbox"/>	Standard seat, XL (SD 42 ... 50 cm, SW 50 ... 56 cm) Add. charge

Settings for the standard seat	
Back angle	<input type="checkbox"/> -9/1/11/21° <input type="checkbox"/> 0/10/20/30°
Not available with electric back angle adjustment!	
Seat inclination SI	<input type="checkbox"/> 0° <input type="checkbox"/> 3° <input type="checkbox"/> 6°
3° and 6° seat inclination not available with seat height 45 cm, electric seat tilt, and seat height adjustment!	
Seat depth SD in cm	<input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50
Seat width SW in cm cm
Lower leg length cm (25 ... 54 cm) (without seat cushion)
Seat height SH in cm	<input type="checkbox"/> 40* <input type="checkbox"/> 45 <input type="checkbox"/> 50 <input type="checkbox"/> 55
*Only with standard seat, small, footrests options EB 01 and EB 07 and without any electrical option possible. (For 3° seat inclination, without seat cushion, dependent on LL, not available with electric seat height adjustment, SH minus LL equals ground clearance, a minimum of 8 cm is recommended)	

Contour seat		
EC 03	<input type="checkbox"/>	Contour seat, size 1 (SD 36 ... 44 cm, SW 38 ... 42 cm) Add. charge
EC 04	<input type="checkbox"/>	Contour seat, size 2 (SD 40 ... 48 cm, SW 43 ... 48 cm) Add. charge
EC 21	<input type="checkbox"/>	Contour seat, XL (SD 40 ... 48 cm, SW 50 ... 56 cm) Add. charge

Settings for the contour seat	
Seat depth SD in cm	<input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48
Seat width SW in cm cm
Lower leg length LL cm (28 ... 54 cm)
Seat height SH in cm	<input type="checkbox"/> approx. 52 <input type="checkbox"/> approx. 57
(For 0° seat inclination, with seat padding, dependent on LL, not available with electric seat height adjustment) (SH minus LL equals ground clearance, a minimum of 8 cm is recommended)	

Recaro seat			
EC 40	<input type="checkbox"/>	Recaro n-Joy (manually adjustable backrest)	Add. charge
EC 41	<input type="checkbox"/>	Recaro n-Joy-Plus (adaptable on site, manually adjustable backrest, seat depth SD 38 ... 56 cm)	Add. charge

Modular Recaro seating system			
EC 42	<input type="checkbox"/>	Backrest LX (manually adjustable)	Add. charge
EC 43	<input type="checkbox"/>	Backrest LT (lateral supports with different setting possibilities, manually adjustable)	Add. charge
EC 44	<input type="checkbox"/>	Seat bottom F (flat padded, seat bottom sewn to lateral guides)	Add. charge
EC 45	<input type="checkbox"/>	Seat bottom X (contoured lateral guides with seat extension)	Add. charge
EC 46	<input type="checkbox"/>	Seat bottom W (flat padded, with seat extension)	Add. charge

Settings for Recaro seat/modular Recaro seating system										
Seat depth SD in cm	<input type="checkbox"/>	38	<input type="checkbox"/>	40	<input type="checkbox"/>	42	<input type="checkbox"/>	44	<input type="checkbox"/>	46
(only for Recaro n-Joy-Plus)	<input type="checkbox"/>	48	<input type="checkbox"/>	50	<input type="checkbox"/>	52	<input type="checkbox"/>	54	<input type="checkbox"/>	56
Lower leg length LL cm (28 ... 54 cm)									
Seat height SH in cm	<input type="checkbox"/>	approx. 52		<input type="checkbox"/>	approx. 57					
<small>(For 0° seat inclination, with seat padding, dependent on LL, not available with electric seat height adjustment) (SH minus LL equals ground clearance, 8 cm is recommended) SH approx. 52 not available with electric seat tilt!</small>										
All Recaro seats include headrest and back bag.										

Chassis options			
EA 13	<input type="checkbox"/>	140 kg/308 lbs chassis (standard equipment)	-
EA 16	<input type="checkbox"/>	200 kg/440 lbs chassis (incl. single-panel aluminium footplate and puncture-proof tyres, reinforced side panel and seat frame)	Add. charge

Chassis accessories			
ES 03	<input type="checkbox"/>	Rear marker plate (required in Germany for road traffic permit)	Add. charge

Seat cushion (only for standard seat)			
EC 06	<input type="checkbox"/>	Seat cushion, black cover (standard equipment)	-
EC 10	<input type="checkbox"/>	Seat cushion with incontinence cover	Add. charge
EC 07	<input type="checkbox"/>	Without cushion	Price reduction
EC 11	<input type="checkbox"/>	Terra seat cushion (up to 140 kg/309 lbs)	Add. charge
EC 12	<input type="checkbox"/>	Terra Aquos seat cushion (up to 140 kg/309 lbs)	Add. charge
EC 13	<input type="checkbox"/>	Aerius seat cushion (up to 140 kg/309 lbs)	Add. charge
EC 14	<input type="checkbox"/>	Aerius Active seat cushion (up to 150 kg/330 lbs)	Add. charge
EC 17	<input type="checkbox"/>	Terra Flair seat cushion (up to 150 kg/330 lbs)	Add. charge
Please refer to our Floam seating system brochure for information on additional seat cushions.			

Seat tilt			
EC 50	<input type="checkbox"/>	Electric seat tilt 30° with centre of gravity shifting (seat inclination 0°, max. load capacity 200 kg/441 lbs)	Add. charge
EC 54	<input type="checkbox"/>	Electric seat tilt 45° with centre of gravity shifting (seat inclination 0°)	Add. charge

Seat height adjustment			
EC 55	<input type="checkbox"/>	Electric seat height adjustment (up to 40 cm, seat inclination 0°, standard equipment on C1000 DS/H) max. load capacity 130 kg/286 lbs!	Add. charge
EC 60	<input type="checkbox"/>	Electric seat height adjustment "Heavy Duty" (up to 30 cm, seat inclination 0°)	Add. charge

Back adjustment		
ED 50	<input type="checkbox"/>	Electric back angle adjustment 30° (not available with Recaro n-Joy and Recaro n-Joy-Plus) Add. charge
(not available in combination with contour seat, small, SD 36 and 38 cm and contour seat, large, SD 40 and 42 cm)		
Back version (only for standard seat)		
ED 01	<input type="checkbox"/>	Back upholstery (standard equipment) -
ED 02	<input type="checkbox"/>	Without back upholstery Price reduction
Belts		
EP 51	<input type="checkbox"/>	Static four-way chest strap kit (only for Recaro seat/seating system) (not available in combination with seat tilt and seat height adjustment) Add. charge
EP 60	<input type="checkbox"/>	Chest belt with hook-and-loop closure, installed at chest height (only for Recaro seat/seating system) Add. charge
EP 65	<input type="checkbox"/>	Lap belt with buckle (standard equipment) -
(Lap belt or four-way chest strap required in Germany for road traffic permit)		
Seating system accessories		
EC 20	<input type="checkbox"/>	Abductor, removable (only for standard seat, not available in combination with electric seat tilt and seat height adjustment) Add. charge
EC 30	<input type="checkbox"/>	Flip-down abductor (only for standard seat, not available in combination with electric seat tilt and seat height adjustment) Add. charge
EE 30	<input type="checkbox"/>	Adapter for headrest installation kit (only for standard and contour seat) Please refer to our accessories catalogue for information about headrests. Add. charge
Footrests		
EB 01	<input type="checkbox"/>	Individual footrests, with plastic footplates (standard equipment) - Available with Standard seat XL and Contour seat XL with max. load capacity of 140 kg/309 lbs; LL up to 50 cm with standard seat.
EB 02	<input type="checkbox"/>	Individual footrests, with aluminium footplates (continuously angle-adjustable) Add. charge Available with Standard seat XL and Contour seat XL with max. load capacity of 140 kg/309 lbs; LL up to 50 cm with standard seat.
EB 03	<input type="checkbox"/>	Single-panel footrest with aluminium footplate Add. charge (continuously angle-adjustable, LL up to 50 cm with standard seat)
EB 04	<input type="checkbox"/>	XXL individual footrests with aluminium footplates Add. charge (continuously angle-adjustable, LL up to 50 cm with standard seat)
EB 07	<input type="checkbox"/>	Without footrest Price reduction
EB 20	<input type="checkbox"/>	Mechanically elevating footrests with length adjustment and aluminium footplates, angle-adjustable, max. load capacity 140 kg/308 lbs! Add. charge (LL from 27 cm for standard seat, from 32 cm for Recaro seat/seating system and contour seat) Cannot be combined with standard seat, XL, and contour seat, XL.
EB 50	<input type="checkbox"/>	Electrically elevating footrest with length adjustment and aluminium footplates Add. charge (up to 140 kg/308 lbs incl. individual aluminium footplates, over 140 kg/308 lbs incl. single-panel aluminium footplate; LL from 27 cm for standard seat, from 32 cm for Recaro seat/seating system and contour seat)
Footrest accessories		
EB 80	<input type="checkbox"/>	Lateral leg pad for electrically and mechanically elevating footrest (pair) Add. charge only for codes EB 20 & EB 50
EB 70	<input type="checkbox"/>	Fixation straps for aluminium footplates (included with mounting accessories) Add. charge

Armrest			
EE 10	<input type="checkbox"/>	Armrest with clothing protector (standard equipment, only for Standard and Contour seats)	-
EE 33	<input type="checkbox"/>	Extended armrest with cut-out (telescoping, only included with Recaro n-joy/njoy+)	Add. charge
EE 40	<input type="checkbox"/>	Extended armrest with clothing protector (telescoping, with cut-out for back tube, only for Standard and Contour seat)	Add. charge

Armrest accessories			
EE 34	<input type="checkbox"/>	Channel forearm pad adapter, fixed, mounting side left	Add. charge
EE 24	<input type="checkbox"/>	Channel forearm pad adapter, fixed, mounting side right	Add. charge
EE 35	<input type="checkbox"/>	Channel forearm pad adapter, swivelling, mounting side left	Add. charge
EE 25	<input type="checkbox"/>	Channel forearm pad adapter, swivelling, mounting side right	Add. charge
EE 36	<input type="checkbox"/>	Arm strap, mounting side left	Add. charge
EE 26	<input type="checkbox"/>	Arm strap, mounting side right	Add. charge

For information on channel forearm pads and hand pads, please refer to our accessories catalogue!

Colour of steering links			
026	<input type="checkbox"/>	Sparkling iron effect light	-
003	<input type="checkbox"/>	Jet black RAL 9005	-
020	<input type="checkbox"/>	Saturn blue metallic	Add. charge
016	<input type="checkbox"/>	Silver metallic	Add. charge
004	<input type="checkbox"/>	Ice blue	Add. charge
018	<input type="checkbox"/>	Signal red RAL 3001	Add. charge
002	<input type="checkbox"/>	Anthracite metallic	Add. charge
025	<input type="checkbox"/>	Stainless steel silver	Add. charge
027	<input type="checkbox"/>	Shimano matte	Add. charge
017	<input type="checkbox"/>	Traffic yellow RAL 1023	Add. charge
013	<input type="checkbox"/>	Jet black RAL 9005 matte	Add. charge
019	<input type="checkbox"/>	White pearl	Add. charge

Tread pattern			
EG 01	<input type="checkbox"/>	Ribs	-
EG 02	<input type="checkbox"/>	Lugs	-
EG 50	<input type="checkbox"/>	Caster wheels with ribs, drive wheels with lugs	-

Tyre colour			
EG 03	<input type="checkbox"/>	Grey	-
EG 04	<input type="checkbox"/>	Black (only with lugs)	-

Tyres			
EG 55	<input type="checkbox"/>	Pneumatic tyres (standard equipment)	-
EG 60	<input type="checkbox"/>	Puncture-proof tyres for all wheels	Add. charge
EF 06	<input type="checkbox"/>	Puncture-proof tyres for caster wheels	Add. charge

Drive unit			
ET 01	<input type="checkbox"/>	6 km/h drive unit (standard equipment)	-
ET 10	<input type="checkbox"/>	10 km/h drive unit (4-pole motors, road traffic permit included, lap belt or four-way chest strap, lighting and rear marker plate required, only available in combination with drive wheel suspension)	Add. charge
ET 06	<input type="checkbox"/>	14 km/h drive unit (not available with 200 kg/441 lbs chassis, road traffic permit included; lap belt or four-way chest strap and rear marker plate required)	Add. charge

Information



EU 40 Calibratable joystick



EU 23 – EU 28 Joystick options



EB 80 Lateral leg pad



EP 51 Four-way chest strap



EB 03 Single-panel footrest



EB 20 Elevating footrest



EB 02 Individual footrests



EC 03/EC 04 Contour seat



EC 40 Recaro N-joy



EC 42 Backrest LX



EC 43 Backrest LT



EC 44 Seat bottom F



EC 45 Seat bottom X



EC 46 Seat bottom W

Batteries			
EU 01	<input type="checkbox"/>	Two 63 Ah (C5) gel batteries (standard equipment)	-
EU 03	<input type="checkbox"/>	Without	Price reduction
EU 05	<input type="checkbox"/>	Two 94 Ah (C5) gel batteries	Add. charge
EU 12	<input type="checkbox"/>	Two 63 Ah (C5) AGM batteries	Price reduction

Battery charger			
EU 91	<input type="checkbox"/>	Battery charger (standard equipment)	-
EU 02	<input type="checkbox"/>	Battery charger without fan (low-noise) for gel batteries (not available with 94 Ah)	Add. charge
EU 92	<input type="checkbox"/>	Without battery charger	Price reduction

Control unit/control panel			
EU 85	<input type="checkbox"/>	Standard control panel (standard equipment)	-
EU 86	<input type="checkbox"/>	Control panel with two external buttons Please indicate installation version: <input type="checkbox"/> Piko button (diameter 50 mm, red/black) <input type="checkbox"/> Piko button (diameter 30 mm, red/black)	Add. charge
EU 87	<input type="checkbox"/>	Mid-tray control (MTC) electrically swivelling (charging/programming box required) Please indicate button type: <input type="checkbox"/> Piko button (diameter 50 mm, red/black) <input type="checkbox"/> Piko button (diameter 30 mm, red/black)	Add. charge

Control panel holder (not for mid-tray control)			
EU 19	<input type="checkbox"/>	Standard control panel holder (standard equipment)	-
EU 20	<input type="checkbox"/>	Control panel holder, swing-away in parallel	Add. charge
EU 18	<input type="checkbox"/>	Control panel holder, swing-away, height-adjustable	Add. charge
EU 29	<input type="checkbox"/>	Control panel holder, downwards adjustable	Add. charge

Control panel holder mounting side			
	<input type="checkbox"/>	Right	
	<input type="checkbox"/>	Left	

Tray-side mounting materials for mid-tray control (MTC)			
	<input type="checkbox"/>	Without	
EU 871	<input type="checkbox"/>	Tray with holder	Add. charge
EU 873	<input type="checkbox"/>	Holder only	Add. charge

External buttons for mid-tray control (MTC)			
	<input type="checkbox"/>	Without	
EU 874	<input type="checkbox"/>	Piko button, diameter 30 mm, red	Add. charge
EU 875	<input type="checkbox"/>	Piko button, diameter 30 mm, yellow	Add. charge
EU 876	<input type="checkbox"/>	Piko button, diameter 30 mm, black	Add. charge
EU 877	<input type="checkbox"/>	Piko button, diameter 50 mm, red	Add. charge
EU 878	<input type="checkbox"/>	Piko button, diameter 50 mm, yellow	Add. charge
EU 879	<input type="checkbox"/>	Piko button, diameter 50 mm, black	Add. charge

One button, "swivel", required; connection of max. 3 buttons, "swivel", "ON/OFF" and "mode" possible.

Charging/programming box for mid-tray control (MTC) (includes 1.5 m BUS cable)			
EU 56	<input type="checkbox"/>	Yes	Add. charge

Activation of the drive-away lock	
<input type="checkbox"/>	Yes
<input type="checkbox"/>	No

Activation of the electronic steering lock	
<input type="checkbox"/>	Yes
<input type="checkbox"/>	No

Special controls		
EU 07	<input type="checkbox"/>	Push-button module (control of max. 10 electric options without LCD) Add. charge
EU 98	<input type="checkbox"/>	Programmable LCD , IR environment control and LED control included Add. charge
EU 97	<input type="checkbox"/>	Programmable LCD , IR environment control, Bluetooth and LED control included Add. charge

The following special controls require an LCD:

EU 09	<input type="checkbox"/>	Chin control with mini joystick, electrically swing-away Add. charge (standard control panel included) Please indicate: <input type="checkbox"/> Soft ball with elastomer <input type="checkbox"/> Soft ball with spring, hard <input type="checkbox"/> Soft ball with spring, soft Please indicate installation version: <input type="checkbox"/> Chin swivel arm, left <input type="checkbox"/> Chin swivel arm, right <input type="checkbox"/> Chin swivel arm overhead, mounting side left (multifunction switch MFS right) <input type="checkbox"/> Chin swivel arm overhead, mounting side right (multifunction switch MFS left)
	<input type="checkbox"/>	Chin control with midi joystick, electrically swing-away Add. charge (standard control panel included) Please indicate installation version: <input type="checkbox"/> Chin swivel arm, left <input type="checkbox"/> Chin swivel arm, right <input type="checkbox"/> Chin swivel arm overhead, mounting side left (satellite switch right) <input type="checkbox"/> Chin swivel arm overhead, mounting side right (satellite switch left)
EU 10	<input type="checkbox"/>	Foot control, left (standard control panel included) Add. charge
EU 10	<input type="checkbox"/>	Foot control, right (standard control panel included) Add. charge
EU 11	<input type="checkbox"/>	Sip and puff control (only with 6 km/h drive unit and Recaro seat/seating system; includes standard control panel) Add. charge
EU 40	<input type="checkbox"/>	Calibratable joystick, standard, includes multifunction button Add. charge (standard control panel included)
EU 45	<input type="checkbox"/>	Calibratable joystick incl. multifunction button for assembly in tray Add. charge (LCD holder with gooseneck for mounting on tray and standard control panel included)

Joystick top version: Please indicate.

	<input type="checkbox"/>	Control knob with elastomer
	<input type="checkbox"/>	Soft ball with elastomer
	<input type="checkbox"/>	Soft ball with spring, hard
	<input type="checkbox"/>	Soft ball with spring, soft
EU 68	<input type="checkbox"/>	Adapter for 1 – 6 buttons Add. charge (required for 1 – 6 button control; select buttons under U 691 to 704)

The following options require an LCD:

EU 61	<input type="checkbox"/>	IR remote control hands-free telephone including mains adapter Add. charge
EU 63	<input type="checkbox"/>	“Dragger 32” mouse click input assistance (included with Satellite environment control and IR receiver mouse emulation in combination with special controls) Add. charge
EU 59	<input type="checkbox"/>	Adapter cable for left/right mouse click Add. charge <input type="checkbox"/> Piko button (diameter 50 mm, red/black) <input type="checkbox"/> Piko button (diameter 30 mm, red/black)

Special control accessories			
EU 69	<input type="checkbox"/>	Easywave wireless module	Add. charge
EU 67	<input type="checkbox"/>	ON/OFF or EMERGENCY STOP with gooseneck for LCD (select buttons under EU 691 to EU 704)	Add. charge
EU 75	<input type="checkbox"/>	ON/OFF button for LCD module, includes MFS adapter cable (select button under EU 691 to EU 704)	Add. charge
		(ON/OFF or EMERGENCY STOP with gooseneck, or ON/OFF button for LED control required)	
Piko button (diameter 50 mm)			
EU 691	<input type="checkbox"/>	Red	Add. charge
EU 692	<input type="checkbox"/>	Yellow	Add. charge
EU 693	<input type="checkbox"/>	Black	Add. charge
EU 697	<input type="checkbox"/>	Silver	Add. charge
Piko button (diameter 30 mm)			
EU 701	<input type="checkbox"/>	Red	Add. charge
EU 702	<input type="checkbox"/>	Yellow	Add. charge
EU 703	<input type="checkbox"/>	Black	Add. charge
EU 704	<input type="checkbox"/>	Silver	Add. charge
LCD holder			
EU 31	<input type="checkbox"/>	Standard LCD holder (standard equipment for LCD)	-
EU 30	<input type="checkbox"/>	LCD holder, swing-away, height- and angle-adjustable	Add. charge
EU 32	<input type="checkbox"/>	LCD holder with gooseneck	Add. charge
EU 34	<input type="checkbox"/>	LCD holder with gooseneck for tray (included)	Add. charge
Control unit accessories			
EU 15	<input type="checkbox"/>	Attendant control	Add. charge
EU 22	<input type="checkbox"/>	Protective bow for standard control panel	Add. charge
EU 23	<input type="checkbox"/>	Fork for tetraplegics, horizontal (adjustable for hand width 70 ... 90 mm)	Add. charge
EU 24	<input type="checkbox"/>	Stick S80 (special joystick handle 80 mm)	Add. charge
EU 25	<input type="checkbox"/>	Soft ball (black cover)	Add. charge
EU 26	<input type="checkbox"/>	Ball top, small (diameter 43 mm)	Add. charge
EU 27	<input type="checkbox"/>	Ball top, large (diameter 48 mm)	Add. charge
EU 28	<input type="checkbox"/>	Joystick top, flexible, including ball top, large	Add. charge
(Extension by 60 mm; control panel holder, downwards adjustable or swing-away, is recommended)			
EU 35	<input type="checkbox"/>	Soft ball, coloured (diameter 30 mm, 4 pieces)	Add. charge
EU 50	<input type="checkbox"/>	Handheld programming device (dealer version, not subject to discount)	Add. charge
EU 55	<input type="checkbox"/>	Programming software with USB adapter (not subject to discount, certificate required)	Add. charge
EU 681	<input type="checkbox"/>	2-fold adapter cable with Piko button (black) (free assignment of 2 electrical functions) Please indicate function: <input type="checkbox"/> Footrest, left <input type="checkbox"/> Footrest, right <input type="checkbox"/> Seat tilt <input type="checkbox"/> Back angle adjustment <input type="checkbox"/> Seat height adjustment <input type="checkbox"/> Light <input type="checkbox"/> Warning flasher <input type="checkbox"/> Direction indicator, left <input type="checkbox"/> Direction indicator, right <input type="checkbox"/> Horn	Add. charge

EU 682	<input type="checkbox"/>	2-fold adapter cable with Piko button (black) (free assignment of 2 electrical functions) Please indicate function: <input type="checkbox"/> Footrest, left <input type="checkbox"/> Footrest, right <input type="checkbox"/> Seat tilt <input type="checkbox"/> Back angle adjustment <input type="checkbox"/> Seat height adjustment <input type="checkbox"/> Light <input type="checkbox"/> Warning flasher <input type="checkbox"/> Direction indicator, left <input type="checkbox"/> Direction indicator, right <input type="checkbox"/> Horn	Add. charge
---------------	--------------------------	--	--------------------

General accessories			
ES 20	<input type="checkbox"/>	Pannier bag, detachable	Add. charge
ES 22	<input type="checkbox"/>	Accessory pouch (e.g. for battery charger, tool kit)	Add. charge
ES 24	<input type="checkbox"/>	Pocket for mobile phone	Add. charge
ES 32	<input type="checkbox"/>	Cover cap for standard seating units	Add. charge
ES 35	<input type="checkbox"/>	Tray, laterally swing-away (includes mounting accessories; cut-out for control panel must be adapted)	Add. charge
EZ 80	<input type="checkbox"/>	Airman pump	Add. charge
ES 37	<input type="checkbox"/>	External kilometre indicator (only in combination with standard control panel and protective bow)	Add. charge
EA 60	<input type="checkbox"/>	Crutch holder	Add. charge
ES 38	<input type="checkbox"/>	Folding rearview mirror, mounting side left	Add. charge
ES 28	<input type="checkbox"/>	Folding rearview mirror, mounting side right (installation on the control panel side is not possible in combination with holder, swing-away in parallel; not available in combination with holder, swing-way in parallel, and tray)	Add. charge
EZ 52	<input type="checkbox"/>	Tool kit	Add. charge
ES 60	<input type="checkbox"/>	Vehicle transportation kit (required for AMF anchor point system) (only available for max. load capacity of 100 kg/220 lbs or less)	Add. charge

Special options			
24817.001	<input type="checkbox"/>	Special colour two-layer coating	Add. charge
24817.000	<input type="checkbox"/>	Special colour single-layer coating	Add. charge
21187.002	<input type="checkbox"/>	Special dimensions for platform for seating shell	Add. charge
40405.002	<input type="checkbox"/>	Back upholstery for power wheelchairs with special length	Add. charge
20465.011	<input type="checkbox"/>	Contour back with standard seat bottom or Contour seat bottom with standard back	Add. charge
36362.001	<input type="checkbox"/>	Avantgarde2 footplates, locked	Add. charge
36362.000	<input type="checkbox"/>	Avantgarde2 footplates, unlocked	Add. charge
25142.001	<input type="checkbox"/>	Amputation legrest for power wheelchairs	Add. charge
22762.003	<input type="checkbox"/>	Push-button module and control panel on the same side	Add. charge
26156.000	<input type="checkbox"/>	Replace satellite switch, by 2 buttons (chin control)	Add. charge
22569.004	<input type="checkbox"/>	Lip joystick (chin control)	Add. charge
23279.000	<input type="checkbox"/>	Mid-tray control according to drawing	Add. charge
19893.001	<input type="checkbox"/>	e-series wireless stop	Add. charge