

## e-support

### Electric auxiliary drive

Maximum load capacity: 125 kg/276 lbs  
(including wheelchair)

**Article number: 492E24=00000\_C**



#### Standard Equipment

##### Servo auxiliary drive for manual wheelchairs operated via the push rings

- Two rear wheels, complete with integrated hub motor and control electronics
- One rechargeable 6.7 Ah nickel-metal hydride battery
- One battery charger
- Two rear wheel adapters with torque support
- Battery installation on the right side
- Standard equipment includes one anti-tipper, installed on the left or right side
- Two product versions:

**Version A for active wheelchair users:** up to 1.5x amplification of arm strength, optimum support, particularly when driving up- or downhill

**Version B for wheelchair users with reduced arm strength:** up to 3 x amplification of arm strength, optimum support on level terrain and gentle inclines

Order no.:  -----	
-------------------------	--

<input type="checkbox"/> Provide estimate only		<input type="checkbox"/> Order	
<input type="checkbox"/> Attachment for custom fabrication enclosed			
Billing address/customer number			
Delivery address			
Sales person		Name	
		Date	
		Signature	

Amplification	
<input type="checkbox"/>	Product version A: 1.5 x amplification
<input type="checkbox"/>	Product version B: 3 x amplification

Adaptation by	
<input type="checkbox"/>	Adaptation by Ottobock
<input type="checkbox"/>	Adaptation by a certified specialist dealer

Adaptation	
<input type="checkbox"/>	Ex works
<input type="checkbox"/>	Second-hand wheelchair (Ottobock)

For use with wheelchair type (please check)		
	<input type="checkbox"/>	Avantgarde T
	<input type="checkbox"/>	Avantgarde VR
	<input type="checkbox"/>	Start M2 Effect
	<input type="checkbox"/>	Start M5 Comfort
	<input type="checkbox"/>	Start M6 Junior
	<input type="checkbox"/>	Motus
	<input type="checkbox"/>	Avantgarde CS
	<input type="checkbox"/>	Avantgarde CV

Rear wheel size		
	<input type="checkbox"/>	Wheel size 22"
	<input type="checkbox"/>	Wheel size 24"

Tyres		
	<input type="checkbox"/>	Pneumatic tyres (standard equipment)
	<input type="checkbox"/>	Puncture-proof

493U24=ST001	<input type="checkbox"/>	Nickel-metal hydride battery	Add. charge
493U24=ST002	<input type="checkbox"/>	Battery charger	Add. charge
493G24=SK012	<input type="checkbox"/>	Spoke protector, transparent (24")	Add. charge
493G24=SK023	<input type="checkbox"/>	Spoke protector, transparent (22")	Add. charge

Note: Puncture-proof tyres are not available with Avantgarde CS/CV.