

## z50 Power add-on drive

Maximum load capacity 125 kg/276 lbs  
(including wheelchair)

**Article number: 492E39=00000\_C**



### Standard Equipment

- 2 rear wheels with integrated hub motor
- VSI control unit
- Removable control panel
- 12.9 Ah (C5) gel batteries
- Battery charger
- Suitable for combination with stair climber
- Adapter plates

#### Note:

When combining a seat depth of 42 cm or less with 24-inch wheels, the battery case is not flush with the seat edge but projects beyond it.

Order no.:	
-----	

<input type="checkbox"/> Provide estimate only		<input type="checkbox"/> Order	
<input type="checkbox"/> Attachment for custom fabrication enclosed			
Billing address/customer number			
Delivery address			
Sales person		Name	
		Date	
		Signature	


Adaptation by	
<input type="checkbox"/>	Ottobock
<input type="checkbox"/>	Certified specialist dealer

Adaptation	
<input type="checkbox"/>	Ex works <span style="float: right;">included</span>
<input type="checkbox"/>	Second-hand wheelchair (Ottobock) <span style="float: right;">upon request</span>

Use with type of wheelchair (please indicate)	
<input type="checkbox"/>	Avantgarde T
<input type="checkbox"/>	Avantgarde VR
<input type="checkbox"/>	Start M2 Effect*
<input type="checkbox"/>	Start M5 Comfort*
<input type="checkbox"/>	Motus
<input type="checkbox"/>	Avantgarde CS
<input type="checkbox"/>	Avantgarde CV
<input type="checkbox"/>	Ventus

\*Note: The diameter of the z50 can differ from the existent manual wheel (wear and tear, make). When changing between manual wheel and z50 this can result in the fact that the function of the knee lever wheel lock is not correct. We therefore recommend adding a matching manual wheel with the same equipment as the z50 to the order.

## Information



ZA 02 Swinging anti-tipper




ZU 40 Exchangeable battery pack 12.9 Ah (C5)




ZU 41 Exchangeable battery pack 17.5 Ah (C5)



ZU 42 For external charging of the batteries



ZG 90 Spoke protector, transparent



ZG 95 Spoke protector with decor

Control panel installation		
	<input type="checkbox"/>	Right
	<input type="checkbox"/>	Left
Wheelbase		
	<input type="checkbox"/>	Active
	<input type="checkbox"/>	Passive
Rear wheel size		
	<input type="checkbox"/>	Wheel size 22" (pneumatic)
	<input type="checkbox"/>	Wheel size 24" (pneumatic, FLEXEL)
Manual wheel		
	<input type="checkbox"/>	Yes (only in combination with Start M2 Effect/Start M5 Comfort) -
Anti-tipper		
ZA 01	<input type="checkbox"/>	Anit-tipper on both sides from wheelchair (standard) -
ZA 02	<input type="checkbox"/>	Swinging anti-tipper <sup>1), 2)</sup> Add. charge
Batteries		
ZU 01	<input type="checkbox"/>	Two gel batteries 12.9 Ah (C5) (standard) -
ZU 02	<input type="checkbox"/>	Two gel batteries 17.5 Ah (C5) Add. charge
Battery options		
ZU 40	<input type="checkbox"/>	Exchangeable battery pack 12.9 Ah (C5) Add. charge
ZU 41	<input type="checkbox"/>	Exchangeable battery pack 17.5 Ah (C5) Add. charge
ZU 42	<input type="checkbox"/>	Charging adapter Add. charge
Tyre options		
	<input type="checkbox"/>	Standard, pneumatic (standard equipment) -
ZG 02	<input type="checkbox"/>	Tyres, puncture-proof (FLEXEL) -
Spoke protectors		
ZG 90	<input type="checkbox"/>	Spoke protector, transparent Add. charge
ZG 95	<input type="checkbox"/>	Spoke protector with decor Add. charge
General accessories		
ZU 20	<input type="checkbox"/>	Control panel holder, swing-away in parallel <sup>2)</sup> Add. charge
ZU 16	<input type="checkbox"/>	Holder for control for attendant <sup>2)</sup> <input type="checkbox"/> right <input type="checkbox"/> left Add. charge
ZU 23	<input type="checkbox"/>	Fork for tetraplegics, horizontal (adjustable for hand width 70 ... 90 mm) Add. charge
ZU 25	<input type="checkbox"/>	Soft ball (black cover) Add. charge
ZU 26	<input type="checkbox"/>	Ball top, small (diameter 43 mm) Add. charge
ZU 27	<input type="checkbox"/>	Ball top, large (diameter 48 mm) Add. charge
ZU 50	<input type="checkbox"/>	Handheld programming device, dealer version (not subject to discount) Add. charge

<sup>1)</sup> With products from other manufacturers upon request

<sup>2)</sup> Start 3 models:

included with active installation

not required with passive installation