



Otto Bock®

Voyager



Engineering excellence is the best way to describe the new Otto Bock Voyager. Designed to fit like a fine suit, the Voyager can be adjusted to a tailor fit... open frame design, horizontally and vertically adjustable rear axle, the back even locks down when folded. When you sit in a Voyager, you become one with the wheelchair. Easy to adjust, lightweight and built to endure the rigors of an active life.

Various adjustment possibilities and a clean, simple design.



red dot
best of the best

The Voyager: Keeping you dynamic.



Voyager in detail

Open frame construction • Folding, angle-adjustable backrest • Continuously adjustable center of gravity • Rear seat height adjustment in 1 cm increments • Offset can be selected: 0/2/4 cm • Continuous lower leg length adjustment • Optional wheel hub brake

Technical data

Voyager	75° and 85° frame offset
Seat width	34-46 cm
Seat depth	36-48 cm
Back height	25-45 cm
Seat height anterior/posterior	43.3-53 / 40-55.5 cm
Weight	from 9.8 kg
Max. load	120 kg
Back angle	15° to 27°
Lower leg length	36-52 cm

Frame colors

Fiery red RAL 3001	Deep black RAL 9005
Green metallic	Sparkle lemon
Jeans metallic	Sparkle light orange
Blue metallic	Sparkle red
Anthracite metallic	Sparkle blue green
Enchante blue	Sparkle Granny Smith
Brilliant silver	Gold glitter

For precise determination of the colors please order our color swatch.