

Seating Systems



- Foot cups
- Fixed and adjustable back and seating systems for improved positioning in wheelchairs

	Order no.: <hr/>
--	---------------------

<input type="checkbox"/> Provide estimate only		<input type="checkbox"/> Order	
Billing address / customer number			
Delivery address			
Sales person		Name	
	 Date Signature	



Backster



"The Back 2"



"The Back 3"



"The Adjustable Seat"

Backster



Wheelchair back with one-hand operation, including pad and cover

Article no.	Matching wheelchairs with a seat width of	Width of the back shell	Back heights of the shell	Contour depth
[] 481D00=SK640	38 – 42 cm	38 cm	45 cm	11 cm
[] 481D00=SK644	42 – 46 cm	42 cm	50 cm	12 cm
[] 481D00=SK648	46 – 50 cm	46 cm	55 cm	13 cm
[] 481D00=SK685	Lateral supports, fixed, small, piece	for Backster size 1 – 3		
[] 481D00=SK686	Lateral supports, fixed, large, piece	for Backster size 2 – 3		
[] 481D00=LK695	Swing-away lateral support, left	for Backster size 1 – 3		
[] 481D00=RK695	Swing-away lateral support, right	for Backster size 1 – 3		
[] 481D00=SK661	Lateral / lumbar pad set	Size 1		
[] 481D00=SK662	Lateral / lumbar pad set	Size 2		
[] 481D00=SK663	Lateral / lumbar pad set	Size 3		
[] 481D00=SK680	Attachment for head rest mounting kit for Backster			
[] 481D00=SV705	Angle set (four-piece) for back rest angle settings of 5° and 10°			

"The Back 2"



Wheelchair back, contoured, with Floam padding

Article no.	Matching wheelchairs with a seat width of	Width of the back shell*	Back heights of the shell	Inclination angle
[] 476D20=AR15-17	38 – 44 cm	33 – 38 cm	40.5 cm	9° to the front and up to 12° to the rear
[] 476D20=AR18-20	46 – 52 cm	40.5 – 45.5 cm	40.5 cm	9° to the front and up to 12° to the rear
[] 476D20=SK010	Back extension			

*different settings in increments of 2.5 cm are possible; not available in combination with angle-adjustable back

"The Back 3"



Wheelchair back, contoured, with Floam padding and lateral supports in tear-drop shape

Article no.	Matching wheelchairs with seat width	Width of the back shell*	Back heights of the shell	Inclination angle
[] 476D30=ARTD15-17	38 – 44 cm	33 – 38 cm	40.5 cm	9° to the front and up to 12° to the rear
[] 476D30=ARTD18-20	46 – 52 cm	40.5 – 45.5 cm	40.5 cm	9° to the front and up to 12° to the rear
[] 476D20=SK010	Back extension			

*different settings in increments of 2.5 cm are possible; not available in combination with angle-adjustable back

"The Back 4"



Wheelchair back, contoured with back padding of Floam cells and lateral supports in tear-drop shape

Article no.	Matching wheelchairs with seat width	Width of the back shell*	Back heights of the shell	Inclination angle
[] 476D40=AR4TD15-17	38 – 44 cm	33 – 38 cm	40.5 cm	9° to the front and up to 12° to the rear
[] 476D40=AR4TD18-20	46 – 52 cm	40.5 – 45.5 cm	40.5 cm	9° to the front and up to 12° to the rear
[] 476D20=SK010	Back extension			

*different settings in increments of 2.5 cm are possible; not available in combination with angle-adjustable back

"The Adjustable Seat"



Firm seat to hang into the wheelchair

Article no.	Matching wheelchairs with seat width	Plate length	Seat depth	Available seat cushion length	Seat height adjustment*	Inclination angle
[] 476C01=DS14-18	40 – 50 cm	44.5 cm	42 – 46 cm	38 – 51 cm	5 up to + 3.5 cm	up to approx. 11°

*Dependent on cross brace design

Foot caps



Fixation sandals with ankle strap, for positioning in the wheelchair, pair

Article no.	Dimensions L x W	Max. load capacity
[] 476B52=SK010	190 x 89 mm	75 kg
[] 476B52=SK015	230 x 114 mm	125 kg
[] 476B52=SK020	250 x 140 mm	180 kg